SAM[®]
NEO

HANDBUCH
MANUAL

Inhalt

Contents

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
*The **Upper Member** of the articulator consists of **Condylar Housings** posteriorly and Upper **Incisal Guidance Table** anteriorly. Together, they simulate the TMJ and upper anterior teeth guidance.*

*The moving **Lower member** of the articulator with **Condylar Elements** posteriorly and the **Incisal Pin** anteriorly, together, they simulate the mandibular condylar movement and lower anterior teeth guidance.*

This, therefore, provides for a very logical and natural laboratory simulation of the patient movements which are very easy to understand and reproduce.

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The red dot in the center of the lower member condylar element corresponds to a hinging axis rotation center of the mandibular condyles. Also known as "Hinge Axis".

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Above the red dotted/dashed line is the axis-orbital plane (AOP)


The upper anterior incisal table surface along with the posterior condylar element centers of rotation are on the same horizontal plane.

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The lower incisal pin in contact with the upper incisal guidance table is used to simulate the relationship of upper and lower incisal edges of the patient's teeth. Thus promoting normal functional patterns.

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

All current articulators and instruments made by SAM[®] are based on the height „P“ (prosthetic height = 121 mm between upper and lower members). The intercondylar distance is 110 mm. A particular feature of the system is the SAM[®] 15 mm distal offset arrangement of the mounting plates in the articulator. This creates a more normal patient simulation of mounted models in the SAM[®] articulator.

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The SAM[®] NEO and all articulators made by SAM[®] are ARCON Instruments. To simulate the patient's authentic TMJ movements, the lower articulator member condylar ball moves in the upper member condylar housing and this movement corresponds to the TMJ anatomical model. This natural anatomical movement is logically simulated and understandable.

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The SAM® NEO is a patient simulator which allows the end-user to reproduce anatomic patient mandibular movements for examination, diagnosis, and therapy.

Aus dem digitalen Artikulator AXIOCOMP
Excerpt from the digital articulator AXIOCOMP

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The articulator axis orbital (Frankfort Horizontal) plane and all horizontal planes of the articulator are parallel with the table top. Centric relation mounted casts are anatomically related to the patient's cranium with enough space for all types of casts, including those with split casts.

The SAM® NEO articulator has the same functional dimensions as all other SAM® articulators (Arcon principle).

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Small kits are delivered in a translucent carrying case, single articulators and batches of 10 are in the special transport cartons.

In addition, articulator mounted casts can be transported in the carrying cases along with the articulators.

It is recommended that all SAM® instrumentation, parts, etc. should be placed in the designated areas in the foam to ensure safe transport and dust-free storage.

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
Support with rubber insert if the articulator is turned upside down

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
Thumb Screw for Inclination of Condylar Housing

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
Underside place for various mounting plates (here MPS magnet mounting plate)

Kondylargehäuse
Condylar Housing

Optional einstellbarer Inzaltisch
Optional Adjustable Incisal Table

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Articulator Lower Member

Transferbogen- Ohrolive- Montagestift
Face Bow Earpiece Mounting Pin

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Schraubmontageplatte
Screw Type Mounting Plate


Magnet mit Magnettof
Magnet with MPS Magnetic Holding Cup

AXIOSPLIT® System
AXIOSPLIT® System

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Draufsicht / *Top View*


*The centric lock is **open** when the spring is compressed.
The position of the bolt can be fixed with the centric lock screw.*

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
Draufsicht / *Top View*

*The centric lock is **closed** when the spring is fully extended as shown.
The position of the bolt can be fixed with the centric lock screw.*

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Caution: Do not touch the end of the spring bolt (as shown here) with your thumb!

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The centric lock is open when the spring is compressed. The position of the bolt can be fixed with the centric lock screw.

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When the locking screw is open, the pressure spring simulates the muscular strength of the patient's jaw. The spring resistance repeatedly moves the condyle ball into the centric position.

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On the other hand, if the centric lock is completely open, sliding movements can be performed.

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
Both Centric Locking Assemblies must be OPEN in order to make Articulator Movements.

AXIOCOMP®

**Der erste Digitale SAM® Artikulator
First Digital SAM® Articulator**

vom Erfinder / by inventor of

AXIOGRAPH®, AXIOTRON® - AXIOQUICK® Recorder

**Unentbehrlich in Didaktik und Lehre!!
Indispensable in Didactics and Teaching**

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


*often called as **Bennett angle** is a Mediotrusion curve and no straight pathway, as in most mean value articulators. The Bennett guide is already given in the condylar housing and corresponds to the most common patient situation in a healthy TMJ. It cannot be set individually.*

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



The Bennett movement is a downward, inward, forward movement or downwards, forward, inward movement all at the same time of the condyle on the mediotrusion (balancing, non-working) side while the laterotrusion (working, rotating) side moves virtually straight lateral for about 0.3 mm.

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An exaggerated mediotrusion movement is shown here. The upper and lower jaws move in the opposite direction. The canine guidance causes disclusion of the posterior teeth when the mandible moves laterally and/or in a protrusive direction.

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To perform a laterotrusion movement, the locking screw of the spring-loaded bolt is loosened.

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By holding the condylar element firmly against the medial wall of the condylar housing and pushing the condylar element downward, inward, and forward the condylar element guiding on the condylar housing pathway will move the lower member towards the midline to perform the mediotrusion movement.

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Average Condylar Pathway Inclination is 35°. The average value relates skull to the Frankfort Horizontal Plane (Axis Orbital Plane). With Complete Dentures an average value of 20° is recommended.

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Setting the of Condylar Pathway Inclination

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Setting the Condylar Pathway Inclination. The adjustment of the inclination is performed by loosening the posterior Thumb screw on the rear of the upper member of the articulator.

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
The inclination of the Condylar Pathway is adjustable from 15° to 65°.

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
In the SAM® NEO Articulator you can adjust the condylar inclination pathway. The condylar pathway curvature is the equivalent of the SAM® white curvature pathway.

The joint path is most often a curved trajectory. This was reflected in the condylar SAM® NEO.

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
The path shown is the real path of the center of rotation of the temporomandibular joint head of the patient.

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To simulate this as faithfully as possible, the condyle of the envelope of condylar pathway elements (joint path curvature) must follow.

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Insert the Tilt Support Rods into the holes on the back of the articulator posts. They support the articulator in tilted positions in an angle of approx. 40°.

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This improves the view of the anterior teeth. Occasionally apply silicone spray, to maintain the O ring (ART 526002).

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The SAM® NEO is like all SAM® articulators which are designed originally for an upside down position as well as upright. Notice how this improves the look at the occlusion while observing and controlling the function and effect of the condylar elements of the articulator.

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Insert the vertical support rod into the hole behind the thumb screw of the incisal table.

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It supports the upper member of the articulator when the articulator is open.

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This improves the view of the occlusion, and makes it easier to work on the model (e.g. wax up).

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Occasionally apply silicone spray to maintain the support rod and the inner O ring (ART 525002).

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The different SAM® Mounting Plate Systems are available for all SAM® articulators and equipment. All SAM® instruments can be converted without difficulty.

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The yellow mounting plate is secured with a Thumb Screw and this provides maximum retention for assembly of models with small amounts of plaster. The mounting plates are intended for single use and for storage and easy reattachment of the mounted models.


The low cost of the mounting plates indicates this is the procedure of choice. Not to mention the saving of the costly time to try and recycle the plates. The positioning of dental cast in the articulator is carried out very precisely by means of two precision alignment pins.

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Benefit:

the patient model remains from beginning to end in the same articulator

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The MPS Magnetic Mounting Plate System is an alternative to the historical screw mounting plate. This Magnetic Plate System allows for the quick removal and replacement of the models in the articulator. The conversion of screw type to magnetic system can be done easily and quickly with the appropriate kits (MSP 101K and 101K MPS).

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The AXIOSPLIT® system is used for trouble-free and accurate exchange of models between different articulators (for example, old and new devices), in practice and laboratory. So the models need only to be sent. The articulators stay in place.

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A Master Check ensures the synchronizing with all SAM® devices. The mounted models can be interchanged without loss of fitting accuracy in all SAM® devices. You will initially need the AXIOSPLIT® Starter Kit with Master Check, norm incisal pin and related matrix holder and an articulator kit for each additional instrument. The exact procedure and more information is described in our AXIOSPLIT-Manual.

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Benefit: Highest precision and your own synchronizing

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Screw Type Mounting Plates Attachment

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The Screw Type Mounting Plate has a round hole behind the threaded nut and an elongated slot anteriorly to assure precise plate positioning on all SAM® instrumentation. set to the precision mounting plate pins. The elongated hole shows to the front.

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
With one hand the plate is held in position and with the second hand thumb and fore finger the thumb screw is tightened to contact. It is further tightened firmly only 1/4 of a turn. tightened with the knurled screw on the other side.

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The mounting plate is fully seated when it sets firmly on the border areas without gap.

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Take note: The mounting plate is screwed firmly to the articulator, the mounted model may not be set easily and quickly removed on and off as like with the MPS system.

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Note: The mounting plate is disposable and should not be used more than once since the difficulty of removing set gypsum and the cleaning of the mounting plate normally will destroy plate and destroy it's accuracy.

MPS-Magnet Montageplatten-Befestigung

MPS-Magnetic Mounting Plates Attachment

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The adhesion plate is set with the rounded edge downwards into the provided hole.

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
The MPS-mounting plate is attached with the round side containing the slot anteriorly onto the positioning pins.

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Take note: The slightly rounded and shiny side shows downwards when it is inserted into the mounting plate (the rounded edge side shows space for plaster).

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To reproduce or build tooth guidances there should to be an incisal table. It must be able to produce the inclination and direction of the tooth guides.

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
To simulate the anterior teeth and canine guidance authentically, the incisal table is always attached to the upper part of the articulator.

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
The incisal table is located as a guiding surface on the top member of the articulator. The flat plane incisal table is used to secure the vertical dimension and to support the incisal pin and can be individualized with plastic via the lateral wings.

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Adjusting the Adjustable Incisal Table (Accessories)

1 Der Inzisaltisch wird nicht ganz bis zum Anschlag in das Artikulator-Oberteil geschoben.
The incisal table is not pushed completely into the articulator top until it stops until the anterior moveable section has been completely raised.

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The incisal pin should be always located in the lower part of the SAM® articulator.

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
The zero mark on the scale in the middle of the incisal pin is normally the starting position for the dental laboratory. The incisal pin can be fixed in the restorations bite height. The incisal pin can be raised to 11 mm or lowered up to 6 mm.

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Due to the articulator construction the desired value must be set twice at the incisal pin while raising the occlusion, so e.g. in an increase of 2.5 mm in the occlusion the adjustment at the incisal pin is +5 mm.

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When lowering of the bite to 1 mm, the incisal pin should be lowered 2 mm.

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Micro-Adjustable Incisal Pin (Accessories)

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The micro adjustable incisal pin can be adjusted accurately and within 0.05 mm range. When lateral knurled screw is opened the pin can be finely adjusted on the upper black scale of incisal continuously.

This can be finely adjusted in the micrometer range.

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*To **Raise**, turn **Clockwise**.
Rotate **Counterclockwise** to **Lower**.*

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Open the lateral knurled screw slightly (about 1/8 turn).

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Rotate the black micro adjustment screw clockwise until the horizontal mark on the scale is in the center of the pin and shows "2" with the vertical mark on the upper black scale "5".

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
The position is fixed with the lateral knurled screw.

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
The axis orbital plane in the patient corresponds in the articulator to the line between the bearing surface of the incisal table to the incisal pin (1) and the centers of the condylar balls (hinge axis) (2).

*Cranium-related mounted models in the articulator are an essential part of diagnosis and therapy.
The correct position of the maxillary model to the reference plane is an absolute prerequisite for model mounting
with a centric registration, but also for setting condylar pathways and even when using mean values.*

*The Frankfort Horizontal (FH) Plane represents a reference plane for the patient and the articulator.
The reference points of the Frankfort horizontal (Porion and orbitale) are identical to the axis-orbital plane (AOE).*

*The scissors mechanism of the SAM[®] anatomic facebow creates and allow for an optimal
positioning of the transfer fork (from a functional and aesthetic point of view).*

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*Please find the exact procedure and the correct technique to attach the facebow and the proper maintenance and care in our AXIOQUICK[®] facebow manual.
This manual can be provided on request by SAM[®] or can be downloaded as PDF file.*

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

The transfer stand SE is used for an easy and fast maxilla mounting in all SAM® articulators and laboratory cast mounting instruments. After the transfer the anatomic facebow can remain in the treatment room, only the transferfork assembly is needed for the model transfer. All versions of the SAM® transfer fork assemblies are compatible with the transfer stand AX. With the transfer stand the upper jaw position of the patient can be quickly and precisely transferred to the articulator.

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A more detailed manual for model mounting can be provided on request by SAM® or can be downloaded as PDF file.

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
Bonwill Triangle and Model Mounting with Rubber Band

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The Bonwill triangle is understood to mean an imaginary triangle in the dentistry, the vertices of which form the lower incisor point (point of contact of the lower central incisors - 24 and 25 in America) and the centers of the two mandibular condyles. This triangle is normally equilateral and its side length is about 10.5 cm.

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This triangle can be simulated in mounting in the articulator according to the average value principle by using a rubber band.

For this purpose, a suitable rubber band is clamped in the ring indentations in the black vertical posts grooves of the NEO articulator and on the zero point of the incisal pin.

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This corresponds roughly to the chewing plane. Meanwhile, the cusp tips can be optically aligned on the imaginary plane which would result from the strained rubber band triangle. The incisal point of the model should be aligned along the tip of the rubber band. It is about 9.5 cm from the back of the triangle.

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Occlusal Plane Indicator (Accessories)

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The occlusal plane indicator set (ART 385) with a plane and a curved plate with Speescher curve can be used for the average-maxillary model mounting or for a prosthetic set up.

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
The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm.

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The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm The angle scale can be adjusted between 0° - 45°, we recommend that the inclination is usually between 0° - 15°.

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


is used in the production of full dentures and as an aid for a setting up of the maxillary teeth. It provides a Spee Curve. The marks on the template indicate auxiliary lines as the middle of jaw ridge. In the loosened state, the cap may also adjust the inclination angle of the remaining teeth can be determined as occlusion plane.

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Cleaning and Maintenance of the Articulator

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The small centric screw on condylar housing, and the set screws on the incisal pin and table should be occasionally cleaned with mild grease solvents and then lubricated with e.g. Multi-use lubricant (ART 193).

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The SAM® NEO is an easy to care for articulator. Plaster residues can be removed with water, soap and a soft brush.

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To remove greasy stains, we recommend mild grease solvents. Condylar ball and housing should always be clean and dry (oil and grease free).

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Accessories and Consumable supplies for SAM® NEO

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Order-Nr. **ATB 334**

Acrylblock
transfer fork support, acrylic XL

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These enduser training take place at the Center for Dental Education in Gauting near Munich. For international guests, we offer courses in English. If you want another language, we can organize an interpreter for extra charge. Please ask us in time. Snacks and drinks are provided

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These enduser trainings can be held in your laboratory, in your surgery or in suitable locations in Germany in German and English. The price includes the training fee, the travel expenses and possibly overnight expenses for the instructors. Course times and dates should be agreed at an early stage.

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These enduser trainings take place in your laboratory, your surgery or in suitable premises in your area within Europe and Switzerland. If you are looking for a specific training, then you will be trained by an "external specialist", a dentist who is familiar with the SAM[®] system and can provide you with case-related explanations in detail.

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We offer enduser training throughout the world to ensure the best possible work for our customers with our products. We are pleased to provide you with an internal employee or an external specialist. The training is held in English, but you can arrange an interpreter for your desired language.

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Are you interested in one of our standard courses? Would you like an individual offer? Contact directly SAM[®] or ask your local distributor.

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SAM® Präzisionstechnik GmbH ist
*SAM® Präzisionstechnik GmbH is