


# MRI MODELL REPOSITIONING INSTRUMENT

This device is used in orthognathic analysis and treatment planning. Three dimensional repositioning of either upper or lower mounted casts is possible with accurate measurements.




## MRI 220K

### MRI Dübelkit

beinhaltet Dübelsetzer (MRI 212), Dübelbohrer (MRI 211) und 200 Dübel (MRI 210)

### MRI plastic sleeve insert kit

includes sleeve insert tool (MRI 212), drill bit (MRI 211), and 200 plastic sleeve inserts (MRI 210)



## MRI 300

### Modell-Repositionierungs-Instrument

nach Prof. Dr. Rainer Schwestka-Polly - für die Operationsplanung in der Kiefer- und Gesichtschirurgie, zur dreidimensionalen Verlagerung des Ober- oder Unterkiefermodells, mit MRI Hülsen-Kit II (MRI 220K)

### model repositioning instrument

according to Prof. Dr. Rainer Schwestka-Polly - used with SAM® 3 articulator mounted casts for diagnosis, analysis and treatment planning in maxillofacial orthopedics, articulator mounted casts can be moved 3D without sectioning, includes MRI sleeves kit II (MRI 220K)



## MRI II Articulator H Kit

nach Prof. Dr. Rainer Schwestka-Polly - für noch mehr Platz bei der Modellmontage - beinhaltet SAM® 3 H (ART 560 / ART 575M) und Modell-Repositionierungs-Instrument (MRI 300)

according to Prof. Dr. Rainer Schwestka-Polly - for even morespace when mounting includes SAM® 3 H (ART 560 / ART 575M) and model repositioning instrument (MRI 300)






MRI 366 MK

## MRI 351K

für Schraubmontageplatten-Systeme / for systems with screw type mounting plates

## MRI 366 MK

für MPS-System / for MPS-system



## MRI II Articulator Kit SAM 3H

beinhaltet / includes

- 01** SAM® 3 Artikulator  
*SAM® 3 articulator*
- 02** Modell-Repositionierungs- Instrument  
(MRI 300)  
*model repositioning instrument (MRI 300)*
- 03** AXIOQUICK® Transferbogenkit - AX  
(ATB 390K)  
*AXIOQUICK® transferbow kit - AX  
(ATB 390K)*
- 04** Teleskop-Bissgabelstütze (ATB 336)  
*telescopic transfer fork support (ATB 336)*
- 05** Magnetsockel für Bissgabelstütze (ATB 338 / ATB339)  
*magnetic block for transfer fork support  
(ATB 338 / ATB 339)*
- 06** Transferstand AX (ATB 398)  
*transfer stand AX (ATB 398)*
- 07** Montagestand (MOH 560) mit 30  
Elastikstäben  
*red mounting stand (MOH 560) with 30  
flexible plastic rods*
- 08** 20 gelbe Montageplatten  
(Schraub- oder Magnetsystem) /  
*20 yellow mounting plates  
(screw or magnetic type)*
- 09** gelber Stapelbox (ART 599)  
*yellow carrying case (ART 599)*





**Rainer Schweska-Polly, DMD**  
Senior Staff Member  
Department of Orthodontics

**Dietmar Kubein-Meesenburg, DMD**  
Professor and Chairman  
Department of Orthodontics

**Hans G. Luhr, MD, DMD**  
Professor and Chairman  
Department of Maxillofacial  
Surgery

Georg-August-University  
Göttingen  
Germany

**Reprint requests:**

Dr Rainer Schweska-Polly  
Department of Orthodontics  
Georg-August-University  
Robert-Koch-Strasse 40  
D-37075 Göttingen  
Germany  
Fax: 49 0551 39 83 50

## Techniques for achieving three-dimensional positioning of the maxilla applied in conjunction with the Göttingen concept



*Three-dimensional repositioning of the maxilla is possible after Le Fort I osteotomy. However, the preoperatively planned and desired position of the maxillary dental arch often cannot be sufficiently achieved during actual surgery, and deviations in the sagittal and vertical dimensions are common. To reduce these errors, the model-repositioning instrument was developed for model surgery in conjunction with the Göttingen concept for orthodontic-surgical treatment with condylar position control. This instrument allows a controlled three-dimensional positioning of jaw segments with three reference points directly on the teeth. The three-dimensional double-splint method combined with a surgical facebow was developed for a controlled three-dimensional positioning of the maxilla during surgery. This instrument and method were applied during treatment of 20 adult patients, and the position of the maxilla before and after surgery was analyzed. It was found that the planned position of the maxillary dental arch could be transferred from model surgery to actual surgery with an accuracy of  $\pm 1$  mm sagittally and vertically. Thus, the application of the Göttingen concept for three-dimensional positioning of the maxilla results in an improvement of accuracy compared with other methods. Furthermore, use of these procedures is easier and less time-consuming during model and actual surgery than are other procedures. (Int J Adult Orthod Orthognath Surg 1998;13:248-258)*

Treatment of patients with severe dentofacial deformities has become routine in modern orthodontics and maxillofacial surgery.<sup>1-4</sup> During treatment the dental arches of the maxilla and mandible are aligned by fixed orthodontic appliances, and the whole dental arch of the maxilla or the mandible can be mobilized and repositioned in the new, preplanned position after Le Fort I osteotomy of the maxilla<sup>5-8</sup> or sagittal split osteotomy of the mandible.<sup>9-12</sup> An increase of bimaxillary surgery has been noted.<sup>13</sup> If

the dentofacial deformity is caused primarily by skeletal malposition of the maxilla to the rest of the skull, such as a maxillary retrognathism or a skeletal open bite, Le Fort I osteotomy is indicated. To bring the condyle into centric relation to the os temporale on each side and to reproduce this position before, during, and after surgery, the Göttingen concept for orthodontic-surgical treatment with condyle position control was introduced about 10 years ago and has since been applied routinely.<sup>14-18</sup>

**Fig 1** (left) Transfer of the maxillary position with an anatomic transfer bow.

**Fig 2** (right) Centric bite registration for mounting the mandibular cast in the articulator.



The model-repositioning instrument and the three-dimensional double-splint method were developed recently for three-dimensional positioning of the dental arches in model surgery and for a three-dimensionally controlled transfer of the position of the maxilla from model surgery to actual surgery, respectively. In this study, the results of the application of these new developments in conjunction with the Göttingen concept are presented.

## Method and materials

### Cast mounting

The casts of the maxilla and mandible are mounted in centric relation in the articulator according to the Göttingen concept. This ensures that the surgical splints can be fabricated with the condyle in centric relation to the os temporale on each side. With use of these splints centric relation can be transferred from the patient before surgery to the patient during surgery. Repositioning of the condyles in centric relation is reproduced during surgery using condylar positioning appliances. This is a necessary condition for a controlled three-dimensional repositioning of jaw segments.



The maxillary cast is mounted by a facebow transfer (Anatomic transfer bow, SAM) in a semiadjustable or fully adjustable articulator (Fig 1). The cast of the mandible is mounted using centric registration (Fig 2). To make the registration as exact as possible and to avoid deformation of the bite registration plate during cast mounting due to discrepan-

cies of the position of the maxillary and mandibular arches, an acrylic resin base is prepared, and registration is taken using four points of Aluwax (Aluwax Dental Products). Using this base for the bite registration, the position of the mandible and both temporomandibular joints can be accurately transferred in centric relation from the patient to the articulator.<sup>19,20</sup>

For model surgery two sets of casts are mounted in two different articulators: According to the Göttingen concept, the presurgical position remains unchanged in the first articulator, and model surgery is performed in the second articulator.<sup>21-28</sup>

### Model surgery with the model-repositioning instrument

Ruler measurements at the base of the articulator above the dental arch and saw cuts along a planned osteotomy line appear to be too inaccurate for repositioning the maxilla during model surgery. Deviations of up to 3.8 mm sagittally and 5.5 mm vertically from the planned position have been observed at the dental arches.<sup>29</sup> These errors occur especially when the maxilla is moved in several dimensions simultaneously. If, for instance, the maxillary dental arch is moved anteriorly and its dorsal part is rotated upward, the incisal edge of a maxillary incisor moves posteriorly away from the planned position. To solve this problem, an appliance for controlled three-dimensional positioning during model surgery was developed.<sup>30</sup> After testing and modification, this appliance,




**Fig 3** (left) Model-repositioning instrument with the maxilla in the preoperative position.

**Fig 4** (right) The maxillary cast is loosened from its base. The tips of the measuring pins are adjusted to the postoperative position and hold the cast.

called the *model-repositioning instrument*, is now in use for routine treatment.

The model-repositioning instrument (SAM) permits simultaneous measurement of the position of three reference points at the dental arch of the cast mounted in the articulator. These reference points are marked on the edge of a maxillary incisor and at the mesiobuccal cusps of the first or second molars bilaterally, and their positions are registered with measuring pins (Fig 3). The measuring pins can be moved in a three-dimensional orthogonal coordinate system. The hinge axis–orbital plane was chosen as the reference plane for this study. This plane is represented by a parallel to the upper surface of the upper part of the articulator. Using these coordinates, the repositioning is exactly defined by the relative sagittal, vertical, and transverse movements of the reference points. Thus, metric repositionings that have been planned clinically or radiographically can be transferred precisely to the cast. For geometric reasons the use of reference points on the teeth leads to a much higher accuracy than reference points on the base of the maxillary cast. Furthermore, the appliance allows an exact transfer of the position of the dental arch from one position to the next. During repositioning the cast is fixed by the tips of the three measuring pins. No osteotomy lines have to be drawn on the base of the cast. In addition, the application of the model-repositioning instrument permits control and adjustment of the position of the maxilla in relation to the mandible.



**Fig 5** The maxillary cast is attached in the desired postoperative position.

After two pairs of casts have been mounted in two different articulators, the second pair of casts is first mounted in the preoperative situation, and the incisal pin and incisal table of the second articulator are removed. Then, the model-repositioning instrument is attached to the upper part of this articulator. The tips of the measuring pins are adjusted to the three reference points directly on the teeth. Thus the coordinates for the reference points in the preoperative position are defined, and they are obtained from measuring scales on the appliance. The maxillary cast is now loosened from its base irrespective of any osteotomy line, and the tips of the measuring pins are subsequently shifted to the coordinates of the desired postoperative situation (Fig 4). The position of the maxillary dental arch in relation to the mandibular dental arch can be controlled and, if necessary, corrected. Subsequently the maxillary cast is attached in the desired new position on the upper part of the articulator (Fig 5). The model-repositioning instrument is removed from the articulator, and the incisal pin and incisal table are reinserted.

Now the maxillary cast is three-dimensionally positioned with respect to the upper part of the articulator, ie, with respect to the chosen reference plane. As a result, the maxillary cast in the second articulator has been moved from the preoperative situation with the condyles in centric relation into the desired postoperative situation while keeping its centric relation.

#### *Splint fabrication according to the three-dimensional double-splint method*

In both articulators self-curing acrylic resin surgical splints are fabricated. According to the Göttingen concept, the first surgical splint is used to transfer and reposition the condyles in the preoperatively obtained centric relation to the patient during surgery. This is important because the preoperative centric relation can differ from the intraoperative centric relation.<sup>31</sup> Without the application of this first splint and condylar position control, centric relation would be incorrect during surgery. The second surgical splint is used to achieve the correct relation of the jaw segments and the occlusion in the planned position.<sup>32-37</sup> In addition, a so-called postoperative splint is fabricated. This splint is constructed like a miniactivator without wires and represents the postoperative dental relation with centric relation of the condyles.

Until now, it has been difficult to transfer the postoperative position of the maxillary dental arch from model surgery to the intraoperative situation with sufficient accuracy. Deviations of the dental arches up to 6.0 mm sagittally and up to 15.0 mm vertically have been observed.<sup>38</sup> These deviations can be explained as follows: After the maxilla is mobilized, the planned dental relationship is obtained using the second surgical splint. However, this dental relationship is only controlled in the sagittal and transverse dimensions. Since the maxillomandibular complex is rotated upward anteriorly and cranially, maintaining centric relation, additional control of the vertical position is necessary. This is usually performed by using one reference point on

the maxilla below the osteotomy line, ie, in the area that is moved during surgery, and a second reference point located either intraorally on the maxilla above the osteotomy line or extraorally on the nose bone.<sup>39-43</sup> Both approaches lead to a displacement of the lower reference point, especially in cases of large anterior advancements or transverse rotations of the maxilla after osteotomy, which makes the vertical placement difficult to control. Deviations for the intraoral approach are higher because the relative error is greater with short measurement distances than with larger extraoral distances.<sup>44,45</sup> For these reasons the exact three-dimensional placement of the maxillary complex in relation to other skeletal structures during surgery was impossible.

To solve this problem, the three-dimensional double-splint method was developed.<sup>46,47</sup> First, the sagittal and transverse placement of the maxillary arch is directed by the surgical splint. Second, the vertical placement is based on the fact that the vertical position of the mandible in relation to a reference plane above the osteotomy line can be exactly reproduced during model surgery and actual surgery. By reproducing a defined mandibular position, the maxilla can therefore be placed three-dimensionally with the aid of surgical splints. Thus, a reference system, which is especially reproducible in the vertical dimension, is chosen for fabrication of these splints in the articulator. This reference system consists of the upper part of the articulator, the condylar boxes, the lower part of the articulator, and the incisal pin.

For preparation of the splints, the bite is opened to the same extent in the preoperative and postoperative articulators, in both cases about 2 or 3 mm. The length of the incisal pin remains constant for the preoperative and postoperative situation (Fig 6). This way it is possible to fabricate a surgical splint for the postoperative situation that represents the planned movements of the maxilla in all three dimensions. The first splint (in the preoperative articulator) is prepared in the articulator representing the preoperative situation. It is applied for positioning




the maxilla and the condyles in centric relation before osteotomy. The second splint (in the postoperative articulator) is fabricated for the controlled three-dimensional movement of the maxilla during actual surgery. For instance, if an upward movement of the maxilla is desired, the second splint will be vertically enlarged in the corresponding way.


According to the Göttingen concept for Le Fort I osteotomies, the two surgical splints are armed with a connection bar for mounting the surgical facebow in the preoperative and postoperative situations on identical positions with respect to the mandibular arch. Before surgery the connection device of the splint in the first articulator, ie, the centric splint, is adapted to the surgical facebow, which is a modification of the system as described by the manufacturer (SAM). The axial pins are then adjusted to centric condylar position of the hinge axis of the articulator (Fig 7).

### *Surgical procedure*

During actual surgery the first splint is inserted, representing the centric condylar position, and temporary intraoperative maxillomandibular fixation is applied. The centric position of each condyle, which is represented by the tips of the axial pins of the surgical facebow, is transferred to the preauricular skin of the patient by a needle puncture with methylene blue (Fig 8). This guarantees that the centric condylar position, which is preoperatively obtained from the patient, will be effectively transferred to the patient during surgery. After the maxilla is mobilized, the second surgical splint is inserted, bearing the connection device for the surgical facebow in an identical position as the first splint. Thus, centric condylar position can be adjusted and controlled. While centric condylar position is maintained, the maxilla is repositioned three-dimensionally, and an osteosynthesis is performed with four minifixation plates and screws (Mini-fixation plate system or Panfix-fixation plate system from LUHR osteosynthesis system, Leibinger). In the case of bimaxillary



**Fig 6** Principle of fabrication of the surgical splints in the articulators according to the three-dimensional double-splint method. The vertical dimension, eg, the length of the incisal pin, is the same in the preoperative and postoperative situations. (Black splint = first splint; striped splint = second splint.)




**Fig 7** Adjustment of the surgical facebow with the preoperative splint in the articulator.



**Fig 8** Application of the surgical facebow during surgery.

**Fig 9** Principle of repositioning the maxilla during surgery according to the three-dimensional double-splint method. The reference length in the calipers is the same in the preoperative and postoperative situation. (Black splint = first splint; striped splint = second splint.)



surgery, the centric position of the condyles during maxillary osteotomy is controlled using a facebow, and the positioning of the condyles during sagittal split osteotomy of the mandible is performed by use of special positioning plates (Mini-fixation plate system from LUHR osteosynthesis system), which are applied bilaterally between the lateral cortex of the mandibular ramus and the anterolateral area of the zygoma.

The three-dimensional placement of the maxilla during surgery is especially based on the fact that the vertical position of the mandible in relation to the skull above the osteotomy line can be exactly reproduced in the preoperative situation as well as in the postoperative situation (Fig 9). The centric condylar position is transferred from the articulator to the patient with the first surgical splint and the facebow. Subsequently a reference point, eg, a drill hole, on the bony structure of the maxilla above the osteotomy line is defined, and its vertical distance to a reference point on the mandible, eg, a premolar bracket, is mea-



sured with calipers. The initially defined distance is arbitrary and depends only on anatomic structures. After the maxillary complex is mobilized, the second surgical splint is inserted, representing the new maxillary position. Temporary maxillomandibular fixation is performed again, and the maxillomandibular complex is rotated upward anteriorly and cranially.

During rotation of the maxillomandibular complex, the desired centric condylar position is controlled and maintained by the facebow, and the vertical dimension is controlled by the calipers, which measure the distance of the two reference points. Since the second splint represents the three-dimensionally controlled position of the maxilla, the distance between the two reference points must be equal to the preoperative distance. This again is controlled with calipers. At this stage, the following points are of main importance: (1) the condyles are in centric relation due to the application of the facebow for positioning control; (2) the position of the mandible in relation to the skull above the osteotomy line is identical to the situation fixed in the first articulator; and (3) the maxilla is three-dimensionally positioned by the second surgical splint.

After osteosynthesis the second surgical splint is removed, and the mandible rotates into the new centric occlusion without any postoperative maxillomandibular fixation. The postoperative splint in combination with light maxillomandibular guiding elastics are inserted.

#### *Cephalometric analysis*

The results of the application of the new developments in the Göttingen concept were analyzed in 20 adult patients. Each patient had received orthodontic-surgical treatment by means of Le Fort I osteotomy alone or in combination with a sagittal split osteotomy. The position of the maxilla for each patient was obtained from lateral radiographs, which were taken preoperatively and postoperatively for routine diagnostic purposes. The sagittal and vertical dimensions of the maxilla were investigated, since




**Fig 10** (left) First surgical splint with three metal balls as reference points (R 11, R 16, and R 27) on the upper surface.

**Fig 11** (right) Lateral radiographs before and after surgery with the first surgical splint in place.

marked deviations from the planned position have been observed in these dimensions. Three reference points in the maxillary dental arch were selected for measurements: the mesial edge of the right central incisor (R 11), and the mesiobuccal cusp of the right first molar (R 16) and the left second molar (R 27). Each was marked with a metal ball (2-mm diameter) inserted in the upper surface of the first surgical splint (Fig 10). Both preoperative and postoperative lateral radiographs were taken with the splint in position in the patient's mouth (Fig 11). In the postoperative situation the splint was lined with wax on the mandibular side.

The computer program WinCeph 4.0 (Compudent) was used to analyze the radiographs. After the tracings of the preoperative and postoperative radiographs were superimposed anatomically with respect to the anterior base of the skull, the achieved displacements of the reference points were measured in the vertical and sagittal dimensions with respect to the hinge axis–infraorbital plane as the reference plane (Fig 12).

The real movement measured on the tracings of the lateral radiograph and the planned movements from model surgery were compared and analyzed using the statistical package STATISTICA (StatSoft). The following values were calculated as essential parameters for the patient groups: (1) the empirical median




**Fig 12** Superimposition of the lateral radiograph tracings and measurement of the movements of reference points R 11, R 16, and R 27. Light lines and teeth = preoperative position; dark lines and teeth = postoperative position. CR = center of rotation; Or = orbitale.


as an estimate for the true median, describing the values for the average patient, and (2) the interquartile difference as a measure for the deviation from the median. Furthermore, the data were presented as Box-Whisker plots. The central value of the box plot represents the median, below which and above which 50% of all values can be observed. The upper and lower boundaries of the box represent the 25% and the 75% quartiles. Finally, the minimal and maximal values are presented.

## Results

The vertical deviation from the planned position of the maxillary complex with reference points R 11, R 16, and R 27 did not exceed 1.5 mm. In the box plot, negative values correspond to



**Fig 13** Deviations of the achieved position of reference points R 11, R 16, and R 27 from the planned position in the vertical dimension. Negative value indicates deviation in cranial direction; positive value indicates deviation in caudal direction.



**Fig 14** Deviations of the achieved position of reference points R 11, R 16, and R 27 from the planned position in the sagittal dimension. Posterior value indicates deviation in cranial direction; anterior value indicates deviation in caudal direction.

a cranial deviation, positive values to a caudal deviation. The median value is zero; hence the errors encountered are due to random effects (Fig 13). In addition to the cranial and caudal deviations, the absolute values of the deviations were calculated as a measure for the absolute error that had occurred. In the vertical dimension, the median for reference point R 11 was 0.3 mm (75% quartile: 0.55 mm; maximum: 1.2 mm); the median for reference point R 16 was 0.25 mm (75% quartile: 0.55 mm; maximum: 1.5 mm); and the median for reference point R 27 was 0.4 mm (75% quartile: 0.8 mm; maximum: 1.2 mm).

For the sagittal dimension, the values for reference points R 11, R 16, and R 27 differed by, at most, 1.7 mm between the achieved position and the planned position. In the box plots, positive values correspond to deviations to the anterior, negative values to posterior deviations. A small systematic tendency of 0.5 mm for the median toward the anterior could be observed (Fig 14), but this tendency is in the range of the accuracy for surgery and for the measurement method. The median of the absolute value of the deviations in the sagittal dimension at R 11 was 0.5 mm (75% quartile: 0.9 mm; maximum: 1.7 mm); the median for R 16 was 0.6 mm; (75% quartile: 0.9 mm; maximum: 1.7 mm); and the median for R 27 was 0.35 mm (75% quartile: 1.3 mm; maximum: 1.2 mm).

Altogether the results show that an orthodontic-surgical treatment of the maxillary complex using the model-repositioning instrument and the three-dimensional double-splint method permits an accuracy of  $\pm 1$  mm in the vertical and sagittal dimensions for the achieved position of the maxilla.

## Discussion

The present study describes the principles and practical application of the model-repositioning instrument and three-dimensional double-splint method for orthodontic-surgical treatment with condylar position control according to the Göttingen concept. As a precondition for the application, the models must be mounted in exact centric relation in the articulator for model surgery, and the centric relation of the condyles must be preserved throughout actual surgery.

The problems with the inaccuracy in measurements at the base of a model during model surgery have been well-known for a long time. Various devices have been developed for a controlled adjustment of single reference points at the maxillary teeth. One method used to precisely position the edges of the maxillary incisors in the sagittal-vertical plane during model surgery is to wind a piece of wire around the incisal pin of the articulator before model surgery. The tip of the wire touches the incisal edge of a maxillary incisor. The wire can be removed

during model surgery and later replaced in its original position. The displacement of the tip of the wire in relation to the position of the maxillary incisal edge demonstrates the movement of the maxillary incisor.<sup>48</sup> Another method is the three-dimensional positioning of the incisal edge of a maxillary incisor by using the model positioning device.<sup>49</sup> Other authors describe methods for successive measurements of different reference points at the teeth by using the modified boley gauge in the articulator<sup>50</sup> or the model surgery platform outside the articulator.<sup>51</sup>

Application of the model-repositioning instrument makes it possible to simultaneously determine the position of three reference points in relation to a reference plane in the articulator and to perform controlled relative movements. In maxillary surgery the dental arch of the mandible defines the position of the maxilla in the sagittal and transverse dimensions. In this case only the vertical dimension needs to be controlled and adjusted by the model-repositioning instrument. In bimaxillary surgery the maxilla is the first mobilized. Since the mandible is mobilized after the maxilla and since its position is also changed, it cannot be used for reference purposes. This increases the difficulties in controlling the position of the maxilla, which now must be positioned with control in all three dimensions. Here, the model-repositioning instrument enables three-dimensional positioning of the maxilla.

The model-repositioning instrument can also be used in surgery with a segmented maxilla. The procedural principles are the same. After the coordinates of the preoperative position are registered, the maxillary cast is loosened from its base and cut in segments. The segments are positioned in the planned postoperative situation corresponding to the dental arch of the mandibular cast, and subsequently the segments are stabilized with plaster. This newly configured maxillary cast is adjusted in the model-repositioning instrument and moved to the postoperative situation. The coordinates for the new position can be read from measuring scales. The application of this new appliance not

only increases accuracy, it is also less time-consuming because the marking of reference lines and cutting along them can be omitted. Furthermore, the maxillary cast can be easily separated from the articulator using a splitcast and moved into the new position.

In actual surgery the precise positioning of the maxilla has often been difficult, especially in regard to its complex three-dimensional repositioning. Problems arose, for instance, when large advancements in the sagittal dimension were planned, when the anterior and posterior regions were moved cranially by different distances, when the midline was rotated, or when the maxillary arch was extended in the transverse direction by an additional paramedian sagittal split in the palatal region. The difficulties of repositioning can be explained by the fact that the positions of the reference points could not be reproduced for the vertical positioning of the maxilla. With use of the three-dimensional double-splint method, a vertical reference line can now be exactly reproduced, since the reference points are above the osteotomy line and at the mandible. This establishes the three-dimensional double-splint method as a new method. The combination of a reproducible reference line, surgical splints, and maintenance of centric condylar position can be seen as the key to successful three-dimensional positioning of the maxilla: In the condylar region, centric relation is obtained with a facebow, while in the anterior region, the defined reference line is reproduced by calipers. If the vertical relation cannot be achieved immediately, it is possible to carefully remove some bone to enable an exact vertical positioning and to create as much bony surface as possible between the osteotomy segments and the rest of the skull. The application of the three-dimensional double-splint method represents a major improvement for the surgeon and introduces additional control when moving the maxilla into the desired position, which in turn means a reduction of the duration of surgery.

The results of the present study show that the model-repositioning instrument

and the three-dimensional double-splint method can be used in combination with standard methods for orthognathic surgery with condylar position control to achieve a controlled and exact positioning of the maxillary complex. This approach is also less time-consuming than other methods. In addition, the stability of the achieved results was found to be high.<sup>52</sup> A necessary precondition for successful application is that the casts, which are mounted by routine for preparation of the osteotomy, must be in centric condylar position. This guarantees that the maxillary segments can be repositioned as planned. The intraoperative positioning of the maxilla with an accuracy of  $\pm 1$  mm in the sagittal and vertical dimensions denotes a major improvement for orthodontic-surgical treatments.

## References

- Bell WH, Proffit WR, White RP. Surgical Correction of Facial Deformities. Philadelphia: Saunders, 1980.
- Epker BN, Fish LC. Dentofacial Deformities. Integrated Orthodontic and Surgical Correction. St Louis: Mosby, 1986.
- Proffit WR, White RP Jr. Surgical-Orthodontic Treatment. St Louis: Mosby, 1991.
- Jaeger A, Schweska-Polly R, Engelke D. Korrektur von ausgeprägten Okklusionsstörungen durch kombiniert kieferorthopaedisch-kieferchirurgische Maßnahmen. In: Kettel W (ed). Deutscher Zahnärztekalendar 1994. Muenchen: Hanser, 1994:147-167.
- Wassmund M. Lehrbuch der praktischen Chirurgie des Mundes und der Kiefer. Leipzig: Barth, 1935.
- Bell WH. Le Fort I osteotomy for correction of maxillary deformities. J Oral Surg 1975;33:412-421.
- Epker BN, Wolford LM. Middle third face osteotomies: Their use in the correction of acquired and developmental dentofacial and craniofacial deformities. J Oral Surg 1975;33:491-508.
- Luhr HG, Kubein-Meesenburg D. Rigid skeletal fixation in maxillary osteotomies: Intraoperative control of condylar position. Clin Plast Surg 1989;16:157-163.
- Trauner R, Obwegeser H. Zur Operationstechnik bei der Progenie und anderen Unterkieferanomalien. In: Trauner R, Obwegeser H (eds). Dtsch Zahn Mund Kieferheilk 1955;23:1-26.
- Dal Pont G. L'osteotomia retromolare per al corruzione della progenia. Minerva Chir 1959;14:1138-1141.
- Hunsuck EE. A modified intraoral sagittal splitting technique for correction of mandibular prognathism. J Oral Surg 1968;26:250-253.
- Luhr HG, Kubein-Meesenburg D, Schweska-Polly R. Bedeutung und Technik der Kiefergelenkpositionierung bei der sagittalen Spaltung des Unterkiefers. Fortschr Kieferorthop 1991;52:66-72.
- Luhr HG, Jaeger A. Indikation, Technik und Ergebnisse der bimaxillaeren Chirurgie. In: Schwenzer N (ed). Die bimaxillaere Chirurgie bei skelettalen Dysgnathien, Sonderband. Stuttgart: Thieme, 1994:1-11.
- Luhr HG. Skelettverlagernde Operationen zur Harmonisierung des Gesichtsprofils—Probleme der stabilen Fixation von Osteotomiesegmenten. In: Pfeifer G (ed). Die Aesthetik von Form und Funktion in der plastischen und Wiederherstellungschirurgie. Berlin: Springer, 1985:87-92.
- Luhr HG, Schauer W, Jaeger A, Kubein-Meesenburg D. Formveraenderung des Unterkiefers durch kieferorthopaedisch-chirurgische Maßnahmen mit stabiler Fixation der Segmente. Fortschr Kieferorthop 1986;47:39-47.
- Kubein D, Luhr HG, Jaeger A, Schauer H-W, von Ehrlich V. Diagnostik der Relation der Kiefergelenke zur Okklusion. Intraoperatives Kontrollverfahren zur Optimierung kieferorthopaedisch-chirurgischer Eingriffe in Verbindung mit der Plattenosteosynthese. Fortschr Kieferorthop 1987;48:267-275.
- Luhr HG. The significance of condylar position using rigid fixation in orthognathic surgery. Clin Plast Surg 1989;16:147-156.
- Luhr HG, Schweska R, Kubein-Meesenburg D. Intraoperative control of condylar position in maxillary osteotomies with rigid skeletal fixation. In: Bell WH (ed). Modern Practice in Orthognathic and Reconstructive Surgery. Philadelphia: Saunders, 1992:628-639.
- Schweska-Polly R, Roese D, Kuhnt D, Walloschek W. Modelloperation und Herstellung von Operationssplinten im Rahmen der gelenkorientierten kieferorthopaedischen Chirurgie mit der "Condylen-Positionierungs-Apparatur"—Part 1: Vorgehen bei Eingriffen am Unterkiefer. Quintessenz Zahntech 1993;19:389-402.
- Schweska-Polly R, Roese D, Kuhnt D, Walloschek W. Modelloperation und Herstellung von Operationssplinten im Rahmen der gelenkorientierten kieferorthopaedischen Chirurgie mit der "Condylen-Positionierungs-Apparatur"—Part 2: Vorgehen bei Eingriffen am Oberkiefer. Quintessenz Zahntech 1993;19:487-499.
- Jung D. Datentransfer zwischen Artikulator und FRS zur Erstellung von diagnostischem Set up und VTO im komplexen Behandlungsfall. Inf Orthod Kieferorthop 1990;22:383-395.
- Ehmer U, Roehling J, Doerr K, Becker R. Calibrated double split cast simulations for orthognathic surgery. Int J Adult Orthod Orthognath Surg 1988;3:223-227.
- Steinhaeuser EW, Janson IM. Kieferorthopaedische Chirurgie—Eine interdisziplinäre Aufgabe. Bd I. Berlin: Quintessenz, 1988.
- Steinhaeuser EW, Rudzki-Janson IM. Kieferorthopaedische Chirurgie—Eine interdisziplinäre Aufgabe. Bd II. Berlin: Quintessenz, 1994.

25. Krenkel C, Lixl G. Ein Modelloperationsgeraet zur Planung und Simulation von OK- und UK-Osteotomien. *Zahnaerztl Prax* 1991;12:471-473.
26. Drescher D, Holtkamp C, Gisbertz M. Die Herstellung von OP-Splints für die orthognathe Chirurgie. *Quintessenz Zahntechnik* 1992;18:507-517.
27. Jost-Brinkmann P-G, Schrinner H-U, Liokumwitsch M. Vorstellung eines Operations-simulationssystems (OSS)—Ein Beitrag zur Verbesserung von Modelloperationen. *Prakt Kieferorthop* 1992;6:89-96.
28. Kaercher H. Three-dimensional craniofacial surgery: Transfer from a three-dimensional model (Endoplan) to clinical surgery—A new technique. *J Craniomaxillofac Surg* 1992;20:125-131.
29. Ellis E III. Accuracy of model surgery: Evaluation of an old technique and introduction of a new one. *J Oral Maxillofac Surg* 1990;48:1161-1167.
30. Schwestka-Polly R, Roeser D, Kuhnt D, Hille K-H. Application of the model-positioning appliance for three-dimensional positioning of the maxilla in cast surgery. *Int J Adult Orthod Orthognath Surg* 1993;8:25-31.
31. McMillen IB. Border movements of the human mandible. *J Prosthet Dent* 1972;27:524-532.
32. Bock D. Die praeoperative Einstellung der Okklusion bei Veraenderungen der Unterkieferlage durch Osteotomie. *Dtsch Zahnaerztl Z* 1968;23:1401-1402.
33. Lindorf HH. Chirurgisch-schaeedelbezugliche Einstellung des Gebisses (Doppelsplintmethode). *Dtsch Zahnaerztl Z* 1977;32:260-261.
34. Ott K. Erweiterte Planung kieferorthopaedisch-chirurgischer Eingriffe. *Oesterr Z Stomatol* 1978;75:332-337.
35. Ellis E. Modified splint design for two-jaw surgery. *J Clin Orthod* 1982;16:619-622.
36. Ripley JF, Steed DL, Flanary CM. A composite surgical splint for dual arch orthognathic surgery. *J Oral Maxillofac Surg* 1982;40:687-688.
37. Somsiri ST. Das Doppelsplintverfahren zur Vorbereitung einer simultanen chirurgischen Lagekorrektur des Ober- und Unterkiefers. *Fortschr Kieferorthop* 1987;48:59-65.
38. McCance AM, Moss JP, James DR. Le Fort I maxillary osteotomy: Is it possible to accurately produce planned pre-operative movements? *Br J Oral Maxillofac Surg* 1992;30:369-376.
39. Turkey TA. Simultaneous mobilization of the maxilla and mandible: Surgical technique and results. *J Oral Maxillofac Surg* 1982;40:96-99.
40. Johnson DG. Intraoperative measurement of maxillary repositioning: An ancillary technique. *Oral Surg Oral Med Oral Pathol* 1985;60:266-268.
41. Nishioka GJ, Van Sickels JE. Modified external reference measurement technique for vertical positioning of the maxilla. *Oral Surg Oral Med Oral Pathol* 1987;64:22-23.
42. Kahnberg K-E, Ridell A. Combined Le Fort I osteotomy and oblique sliding osteotomy of the mandibular rami. *J Craniomaxillofac Surg* 1988;16:151-156.
43. Kahnberg K-E, Sunzel B, Astrand P. Planning and control of vertical dimension in Le Fort I osteotomies. *J Craniomaxillofac Surg* 1990;18:267-270.
44. Van Sickels JE, Larsen AJ, Triplett RG. Predictability of maxillary surgery: A comparison of internal and external reference marks. *Oral Surg Oral Med Oral Pathol* 1986;61:542-545.
45. Stanchina R, Ellis E III, Gallo WJ, Fonseca RJ. A comparison of two measures for repositioning the maxilla during orthognathic surgery. *Int J Adult Orthod Orthognath Surg* 1988;3:149-154.
46. Schwestka R, Engelke D, Kubein-Meesenburg D. Condylar position control during maxillary surgery: The condylar positioning appliance and three-dimensional double splint method. *Int J Adult Orthod Orthognath Surg* 1990;5:161-165.
47. Schwestka-Polly R, Engelke D, Kubein-Meesenburg D. Application of the condylar positioning appliance in mandibular sagittal split osteotomies with rigid skeletal fixation. *Int J Adult Orthod Orthognath Surg* 1992;7:15-21.
48. Ellis E III, Gallo WJ. A method to accurately predict the position of the maxillary incisor in two-jaw surgery. *J Oral Maxillofac Surg* 1984;42:402-404.
49. Schwestka R, Engelke D, Zimmer B, Kubein-Meesenburg D. Positioning control of the upper incisors in orthognathic surgery. Preoperative planning with the "Model Positioning Device" and intraoperative application of the "Sandwich Splint." *Eur J Orthod* 1991;13:367-371.
50. Henry CH. Modified boley gauge for use as a reference plane in orthognathic surgery. *J Oral Maxillofac Surg* 1990;48:535-539.
51. Erickson KL, Bell WH, Goldsmith DH. Analytical model surgery. In: Bell WH (ed). *Modern Practice in Orthognathic and Reconstructive Surgery*. Philadelphia: Saunders, 1992: 154-205.
52. Jaeger A, Redeker K, Merten H-A. Ergebnisse nach kieferorthopaedisch-chirurgischer Korrektur der Angle-Klasse III mit Hilfe isolierter Le Fort I-Osteotomie oder bimaxillaerer Umstellungsosteotomien. Abstractband Jahrestagung der Deutschen Gesellschaft für Kieferorthopaedie, Bremen, 1996.





Rainer Schweska-Polly, DMD  
Senior Staff Member  
Department of Orthodontics  
Georg-August-University  
Robert-Koch-Strasse 40  
D-3400 Goettingen  
Germany

## Significance of the contour of the lateral surface of the maxilla for planning osteotomy lines in orthognathic surgery

*The Le Fort I osteotomy is one of the most frequently used procedures in orthognathic surgery. The preoperative planning of a Le Fort I osteotomy usually involves a two-dimensional lateral radiograph. On the lateral radiograph, the movement of the maxilla in the sagittal and vertical planes can be simulated, and the position of the osteotomy line in relation to the occlusion can be determined. However, the lateral radiograph consists of two dimensions, whereas, during surgery, three dimensions have to be considered. The contour of the lateral surface of the maxilla and its individual regional variations are not reflected on the lateral radiograph. Therefore, this third dimension's influence on treatment planning is investigated. Discrepancies of the position of planned osteotomy cuts are described after the transfer from two to three dimensions. Discrepancies of the positions of reference points and lines are up to 3.3 mm. Therefore, it is not possible to transfer reference points and reference lines from the lateral radiograph to the maxilla nor to move the mobilized maxilla along the planned osteotomy lines in a precise manner. It is possible to indicate the general direction of the osteotomy lines. Results of this study show the necessity for the precise control of three-dimensional positioning of the maxilla during treatment planning, cast surgery, and actual surgery. (Int J Adult Orthod Orthognath Surg 1993;8:191-201.)*

### Introduction

Modern orthodontics and oral and maxillofacial surgery enable a combined approach for correction of dentofacial deformity. The dental arches of the maxilla and mandible are aligned by fixed orthodontic appliances. The maxilla and/or mandible can be mobilized with surgery, adjusted (three-dimensionally) in a new position, and fixed with miniplates and screws. Following surgery, a fine adjustment of occlusion and muscular rehabilitation are completed by orthodontics.

In orthognathic surgery the Le Fort I osteotomy and subsequent three-dimensional positioning of the maxilla is frequently performed.<sup>1-4</sup> This method might be isolated or combined with mandibular osteotomy. Surgery is planned after clinical evaluation and record analysis. On the lateral radiograph, the surgical movement of the maxilla in the sagittal and vertical planes and the position of the osteotomy cuts in relation to occlusion can be planned.

Reference lines are drawn perpendicular to the occlusal plane from the




Fig 1 (left) Determination of the vertical reference lines and the position of the osteotomy cuts on the tracing from the lateral radiograph.




Fig 2 (right) Simulation of the maxillary movement along the osteotomy cuts.



Fig 3 (left) Two-dimensional planning on the lateral radiograph.




Fig 4 (right) Three-dimensional transfer on a skull.


tips of the canine and the mesiobuccal cusp of the first molar to the maxilla. By linking the most cranial points of these reference lines at the level of the nasal floor, the position of the osteotomy cuts can be determined (Fig 1). According to these reference lines, the movement of the maxilla can be simulated in the sagittal and vertical dimensions (Fig 2).

During surgery the reference lines from the lateral radiograph are transferred to the lateral surface of the maxilla using a caliper. Treatment planning usually consists of two dimensions (Fig 3), whereas, during actual surgery, three dimensions must be considered. The contour of the lateral surface of the maxilla and its individual variations are

Fig 5 (left) Lateral radiograph of a studied skull.

Fig 6 (right, above) Individual impression tray for the maxilla.

Fig 7 (right, below) Cast of a maxilla from a studied skull.





not considered on the lateral radiograph, but this contour is important during surgery (Fig 4). This problem is well known.<sup>5,6</sup> However, no study examines the extent of discrepancies when a planned osteotomy line is transferred from two into three dimensions. The aim of this study is to evaluate the discrepancies that arise when an osteotomy line is transferred from a lateral radiograph onto the maxilla using a caliper.

#### Method and materials

For this study, lateral radiographs were taken of 20 skulls (Fig 5). An individual impression tray was fabricated from self-curing resin (Fig 6). Impressions were taken of the maxilla and maxillary teeth with this tray. Casts were fabricated (Fig 7), and cuts were made in the transverse vertical plane at the level of the tips of the canines and the mesiobuccal cusps of the first molars bilaterally (Fig 8). Contour lines of the lateral surfaces of the cuts in the areas



of the molars (Fig 9) and canines (Fig 10) were drawn with a pencil on graph paper. The cuts were so orientated that a line linking the tips of the left and right tooth was parallel to a horizontal line on the graph paper. Reference lines were drawn perpendicular to the occlusal plane from the tips of the canine and mesiobuccal cusp of the first molar to the maxilla. A desired height represented the length of reference line  $S_1$  with end point  $P_1$ . The lengths of  $S_1$  were chosen as tooth length (from the lateral radiograph) plus 2 mm and tooth length plus 4 mm. Therefore, an adequate distance from the tooth's roots was guaranteed during osteotomy, and different positions of osteotomy cuts could be examined (Fig 11).

The reference line  $S_1$  was transferred from the lateral radiograph to the contour lines using a caliper. On the contour lines the base point of  $S_1$  was the tip of the reference tooth; the vertical direction of  $S_1$  ended cranially at point  $P_1$ .  $S_1$  was drawn a second time so that



*Fig 8 (left)* Cuts on the casts in the transverse vertical plane at the level of the (a) mesiobuccal cusps of the first molars and the (b) tips of the canines bilaterally.

*Fig 9 (right)* Contour curves of the lateral surfaces of the cuts in the area of the molars.



*Fig 10 (left)* Contour curves of the lateral surfaces of the cuts in the area of the canines.

*Fig 11 (right)* Reference line  $S_1$  and end point  $P_1$  on a tracing of the lateral radiograph.




Fig 12 Transfer of the planned reference lines and end points from the lateral radiograph to the contour lines.




Fig 13 Transfer of the planned reference lines and end points for measurement and calculation of the discrepancy  $S_2$  minus  $S_1$ , corresponding the lateral radiograph on the contour lines.

its base point ran through the tip of the tooth again, but the end point  $\bar{P}_2$  overlapped the contour line, which represented the lateral contour of the maxilla.

After  $S_1$  was drawn on the contour line,  $S_1$  and  $\bar{P}_2$  were projected onto the sagittal vertical plane. The desired heights of  $S_1$ ,  $\bar{P}_1$ , and  $P_1$  could not be obtained (for geometric reasons) using a caliper due to the inclination of the lateral contour of the maxilla; only the height of point  $\bar{P}_2$ , which lies more caudal than  $\bar{P}_1$  and  $P_1$ , could be determined. (Fig 12). From the base point of  $S_1$  (the tip of the examined tooth) and  $\bar{P}_2$ , a triangle was constructed with the base of  $S_1$  and the sides of  $S_2$  and  $S_3$  (Fig 13). The length of  $S_3$  was measured, and the length of  $S_2$  was calculated using the Pythagorean theorem ( $S_2^2 + S_3^2 = S_1^2$ ). The difference of the lengths of the lines  $S_2$  minus  $S_1$  characterizes the discrepancy between the height of the points  $\bar{P}_1$  or  $P_1$  on the lateral radiograph and the height of point  $\bar{P}_2$  on the contour line, in projection onto the sagittal vertical

plane. The difference of  $S_2$  minus  $S_1$  is negative because  $\bar{P}_2$  is more caudal than  $\bar{P}_1$  or  $P_1$ . The discrepancies are noted in statistical evaluations (mean, standard deviation, minimum, maximum).

When  $S_1$  is drawn on the lateral contour of the maxilla, the desired height of point  $P_1$  (with a distance  $S_1$  from the base point) cannot be achieved. Only the height of point  $\bar{P}_2$  can be determined (with a distance  $S_2$  from the base point) when projected onto the sagittal vertical plane. Therefore, the length of  $\bar{S}_2$  was calculated to obtain the desired height of point  $P_1$  (with a distance  $S_1$  from the base point) when projected onto the sagittal vertical plane—to obtain the exact height from the lateral radiograph. From the base point on the tip of the examined tooth and from the points  $\bar{P}_1$  and  $P_1$ , a triangle was constructed with the base  $\bar{S}_2$  and the sides  $S_1$  and  $\bar{S}_3$  (Fig 14). The length of  $\bar{S}_3$  was measured, and the length of  $\bar{S}_2$  was calculated using the Pythagorean theorem




Fig 14 Transfer of the planned reference lines and end points for measurement and calculation of the discrepancy  $\bar{S}_2$  minus  $S_1$ , corresponding the lateral radiograph on the contour lines.




Fig 15 Statistical evaluations of the discrepancy  $S_2$  minus  $S_1$  for the first molar.

$(S_1^2 + \bar{S}_3^2 = \bar{S}_2^2)$ . The difference of  $\bar{S}_2$  minus  $S_1$  characterizes the value of the length, which must be added to obtain the desired height of  $S_1$  on the skull. The difference of  $\bar{S}_2$  minus  $S_1$  is positive. The discrepancies are noted in statistical evaluations (mean, standard deviation, minimum, maximum).

**Results**

The discrepancy of  $S_2$  minus  $S_1$  for the first molar averaged 0.3 to 1.3 mm, with

a minimum of 0 mm and a maximum of 3.3 mm. The statistical evaluation shows severe discrepancies, which depend on the contour of the maxilla and the length of the examined reference lines—the longer the reference lines, the larger the differences. Differences were also noted between the right and left skull halves. The distribution of the varying values on the right and left sides is random (Fig 15).

The discrepancy of  $S_2$  minus  $S_1$  for the canine averaged 0.2 to 0.3 mm with a minimum of 0 mm and a maximum of

Fig 16 Statistical evaluations of the discrepancy  $\bar{S}_2$  minus  $S_1$  for the canine.





Fig 17 Statistical evaluations of the discrepancy  $\bar{S}_2$  minus  $S_1$  for the first molar.



0.7 mm. Therefore, discrepancies are distinctly less than those for the first molar area (Fig 16).

The discrepancy of  $\bar{S}_2$  minus  $S_1$  for the first molar, depending on the tooth's length, averaged 0.6 to 2.1 mm with a minimum of 0 mm and a maximum of 6.5 mm. The statistical values show that there are severe discrepancies (Fig 17).

The discrepancy of  $\bar{S}_2$  minus  $S_1$  for the canine averaged 0.2 to 0.3 mm with a minimum of 0 mm and a maximum of 0.7 mm (Fig 18).

### Discussion

Because of the anatomy of the lateral contour of the maxilla, reference lines




Fig 18 Statistical evaluations of the discrepancy  $\bar{S}_2$  minus  $S_1$  for the canine.

for osteotomy cuts determined by the lateral radiograph could not be transferred onto the maxilla in a precise manner. It is not possible to exactly move the mobilized maxilla along the planned osteotomy cuts. In the region of the first molar, the discrepancies may be several millimeters and cannot be estimated when the exact anatomy of the patient is unknown. The longer the reference lines are, the bigger the deviations—due to the inclination of the maxilla in the region of the zygomatic buttress. The discrepancies for the first molar and the canine are different. In the region of the canines, the discrepancy is relatively small and can be ignored for clinical purposes. The discrepancies that may exist can change the location of the planned osteotomy cuts, particularly when longer reference lines have to be transferred in the region of the molars. Distinct differences may also exist between the right and left sides. The distribution of the right and left sides may be random, but this only raises inaccuracy on the whole.

In this study discrepancies between reference positions for osteotomy cuts on the lateral radiograph and the maxilla are described quantitatively for the

first time. These discrepancies must be accommodated during treatment planning, cast surgery, and actual surgery. Different approaches have been introduced.

One approach is to transfer reference lines more exactly from the lateral radiograph to the maxilla during cast surgery and actual surgery with the maxillary measuring appliance.<sup>7</sup> This device consists of a flat horizontal measuring table (30 × 10 mm) and a pillar (35 mm) made of metal—the pillar being vertical to the measuring table. The pillar is fastened to the longer side of the measuring table and can be moved along this side. During cast surgery the measuring table is fastened to the occlusal surface of the mounted maxillary cast extending from the canine to the first molar with self-curing resin. The vertical pillar enables the transfer (with a pencil) of vertical and horizontal reference lines from the lateral radiograph to the maxillary cast. The geometric error from a caliper is avoided. During actual surgery the measuring table is positioned on the dental arch of the patient. With the vertical pillar, vertical and horizontal reference lines can be transferred (with a



bur) to the lateral walls of the maxilla. The application of the device appears relatively complicated, because two devices are necessary for the right and the left sides.

The second approach is to make the treatment planning not two dimensional, but three dimensional.<sup>8-13</sup> One procedure utilizes a radiograph in a posteroanterior direction, providing a second plane for treatment planning to obtain information about the anatomy of the lateral contour of the maxilla.

Another procedure utilizes computer tomography. From tomographic scans the anatomy of the skull can be visualized three dimensionally. This data can be integrated to produce plastic models in a milling machine. Using these models, orthognathic surgery can be simulated. This method requires equipment and time. This procedure may be indicated for special treatment planning of severe dentofacial deformities.

In addition to the described discrepancies due to the maxillary contour, the width of the saw cut for the osteotomy and three-dimensional movement of the mobilized maxilla must be considered. This is especially relevant with problems of asymmetry when complex three-dimensional movements are performed, and telescoping appears between the mobilized maxilla and the rest of the skull.<sup>14</sup> The goal is to make a controlled three-dimensional repositioning of the maxillary teeth and the entire maxilla during treatment planning, cast surgery, and actual surgery.<sup>14-28</sup> It is also important to consider and to position the cant of the occlusal plane with special surgical techniques.<sup>29-31</sup> Additional procedures are described in previous literature.<sup>32-60</sup>

The third approach assumes that only the general direction of an osteotomy line is important, but not its exact position. The position of the osteotomy line can run parallel to the occlusal plane or at an angle to it.

Through advanced technologies, the planned position of the maxilla can be exactly transferred from cast surgery to actual surgery. The application of advanced positioning techniques during cast surgery, such as with the model po-

sitioning appliance, and the application of advanced surgical splint technologies, such as the sandwich splint or the three-dimensional double splint in combination with the condylar positioning appliance, represent a closed concept and allow three-dimensional controlled positioning of the maxilla after Le Fort I osteotomies. This concept is independent of an exact positioning of osteotomy cuts, but it requires an exact control in the vertical dimension. This third approach is relatively easy to handle in daily routine.

The anatomy is the basis for the development and critical discussion of clinical procedures. Knowing the exact anatomy of the maxilla will help to avoid errors during planning and treatment of dentofacial deformities. Measurements obtained at the planning stage must be clinically modified if accurate repositioning of the maxilla is critical to a successful outcome.

## References

1. Le Fort R. Etude experimentale sur les fractures de la machoire superieur. *Rev Chir* 1901;23:208-227, 360-379, 479-507.
2. Wassmund M. *Lehrbuch der praktischen Chirurgie des Mundes und der Kiefer*. Leipzig: Meuser, 1935.
3. Bell WH. Le Fort I osteotomy for correction of maxillary deformities. *J Oral Surg* 1975; 33:412-426.
4. Epker BN, Wolford LM. Middle third facial osteotomies: their use in the correction of acquired and developmental deformities and cranio-facial deformities. *J Oral Surg* 1975; 33:491-514.
5. Austermann KH, Bollmann F. Eine Methode zur Bestimmung der guenstigsten Osteotomieebene bei Le-Fort-I-Osteotomien. *Fortschr Kiefer Gesichtschir* 1981;26:121-123.
6. Muehling J, Reuther J. Versorgung von Mittelgesichtstruemerfrakturen mit dem Wuerzburger Titan-Miniplatten-System. *Fortschr Kiefer Gesichtschir* 1991;36:70-73.
7. Wylie GA, Epker BN, Mossop JS. A technique to improve the accuracy of total maxillary surgery. *Int J Adult Orthod Orthognath Surg* 1988;3:143-147.
8. Marsh JL, Vannier MW. The "third" dimension in craniofacial surgery. *Plast Reconstr Surg* 1983;71:759-767.
9. Hirschfelder U, Hirschfelder H. 3-D-Rekonstruktion zur Beurteilung der Morphologie kraniofazialer Strukturen. *Dtsch Zahnärztl Z* 1989;44:187-190.
10. Lambrecht JT, Brix F. Planning orthognathic surgery with three-dimensional models. *Int J Adult Orthod Orthognath Surg* 1990;5:141-144.
11. Marsh JL, Vannier MW. The role of 3-D reformation from CT scans in craniofacial

- surgery. In: Pfeifer G (ed). *Craniofacial Abnormalities and Clefts of the Lip, Alveolus and Palate*. 4th Hamburg International Symposium. Stuttgart: Thieme, 1991:46–48.
12. Ono I, Ohura T, Narumi E, Kawashima K, Matsuno I, Nakamura S, et al. Three-dimensional analysis of craniofacial bones using three-dimensional computer tomography. *J Craniomaxillofac Surg* 1990;20:49–60.
  13. Stoker NG, Mankovich NJ, Valentino D. Stereolithographic models for surgical planning: preliminary report. *J Oral Maxillofac Surg* 1992;50:466–471.
  14. Turvey TA, White Jr RP. Maxillary Surgery. In: Proffit WR, White Jr RP (eds). *Surgical-Orthodontic Treatment*. St Louis: Mosby, 1991:248–263.
  15. Haerle F. Le Fort I osteotomy (using miniplates) for correction of the long face. *J Oral Surg* 1980;38:427–432.
  16. Zuercher A, Hardt N, Steinhäuser EW. Interokklusaler Splint als Rezidivprophylaxe bei totalen Unterkieferosteotomien. *Fortschr Kiefer Gesichtschir* 1981;26:81–82.
  17. Ellis E. Modified splint design for two-jaw surgery. *J Clin Orthod* 1982;16:619–622.
  18. Ripley JF, Steed DL, Flanary CM. A composite splint for dual arch orthognathic surgery. *J Oral Maxillofac Surg* 1982;40:687–688.
  19. Turvey TA, Hall DJ, Fish LC, Epker BN. Surgical-orthodontic treatment planning for simultaneous mobilization of the maxilla and mandible in correction of dentofacial deformities. *Oral Surg Oral Med Oral Pathol* 1982;54:491–498.
  20. Lühr HG. Miniplate fixation of Le Fort I osteotomies [Discussion]. *Plast Reconstr Surg* 1986;78:755.
  21. Somsiri ST. Das Doppelsplintverfahren zur Vorbereitung einer simultanen chirurgischen Lagekorrektur des Ober- und Unterkiefers. *Fortschr Kieferorthop* 1987;48:59–66.
  22. Bell WH, Mannai C, Lühr HG. Art and science of the Le Fort I downfracture. *Int J Adult Orthod and Orthognath Surg* 1988;3:23–52.
  23. Lindorf HH. Osteosynthese durch Schrauben und Miniplatten bei kieferorthopaedischen Operationen. *Inf Orthod Kieferorthop* 1988;20:329–350.
  24. Lühr HG, Kubein-Meesenburg D. Rigid skeletal fixation in maxillary osteotomies. Intraoperative control of condylar position. *Clin Plast Surg* 1989;16:157–163.
  25. Paulus GW. Moderne funktionelle und aesthetische Aspekte bei kieferorthopaedischen Operationen. *Inf Orthod Kieferorthop* 1990;22:33–55.
  26. Wangerin K. Einzeitige bimaxillaere Korrektur extremer Fehlbisse—Vorbehandlung, planung und Operationsmethode mit funktionsstabiler Fixierung im Ober- und Unterkiefer. *Dtsch Z Mund Kiefer Gesichtschir* 1990;14:424–432.
  27. Lühr HG, Schweska R, Kubein-Meesenburg D. Intraoperative control of condylar position in maxillary osteotomies with rigid skeletal fixation. In: Bell WH (ed). *Modern Practice in Orthognathic and Reconstructive Surgery*. Philadelphia: Saunders, 1992:628–639.
  28. Merten H-A, Halling F. A new condylar positioning technique in orthognathic surgery. *J Craniomaxillofac Surg* 1992;20:310–312.
  29. Reyneke JP, Masureik CJ. Treatment of maxillary deficiency by a Le Fort I downsliding technique. *J Oral Maxillofac Surg* 1985;43:914–916.
  30. Munro IR. The self-retained Le Fort I osteotomy. In: Marchac D (ed). *Craniofacial Surgery*. Berlin: Springer, 1987:326–329.
  31. Reyneke JP, Evans WG. Surgical manipulation of the occlusal plane. *Int J Adult Orthod Orthognath Surg* 1990;5:99–110.
  32. Jung D. Datentransfer zwischen Artikulator und FRS zur Erstellung von diagnostischem Set up und VTO im komplexen Behandlungsfall. *Inf Orthod Kieferorthop* 1988;20:383–395.
  33. Bell WH, Proffit WR, White RP. *Surgical Correction of Facial Deformities*. Philadelphia: Saunders, 1980:234–441.
  34. Epker BN, Fish LC. *Dentofacial Deformities: Integrated Orthodontic and Surgical Correction*. St Louis: Mosby, 1986.
  35. Ehmer U, Roehling J, Doerr K, Becker R. Calibrated double split cast simulations for orthognathic surgery. *Int J Adult Orthod Orthognath Surg* 1988;3:223–227.
  36. Steinhäuser EW, Janson I. *Kieferorthopaedische Chirurgie—Eine interdisziplinäre Aufgabe*. Berlin: Quintessenz, 1988.
  37. Krenkel C, Lixl G. Ein Modelloperationsgeraet zur Planung und Simulation von OK- und UK-Osteotomien. *Zahnaerztl Prax* 1991;42:471–473.
  38. Proffit WR, White Jr RP. *Surgical-Orthodontic Treatment*. St Louis: Mosby, 1991.
  39. Jost-Brinkmann P-G, Schrinner H-U, Liukowitsch M. Vorstellung eines Operationssimulationssystem (OSS)—Ein Beitrag zur Verbesserung von Modelloperationen. *Prakt Kieferorthop* 1992;6:89–96.
  40. Kaercher H. Three-dimensional craniofacial surgery: transfer from a three-dimensional model (Endoplan) to clinical surgery—a new technique (Graz). *J Craniomaxillofac Surg* 1992;20:125–131.
  41. Morris JH, Albright JE, Brooks JH, Naito RL. Prediction of ramus fragment relationship in the surgical correction of facial asymmetries. *Int J Adult Orthod Orthognath Surg* 1992;7:245–250.
  42. Ellis III E. Accuracy of model surgery: Evaluation of an old technique and introduction of a new one. *J Oral Maxillofac Surg* 1990;48:1161–1167.
  43. Ellis III E, Gallo WJ. A method to accurately predict the position of the maxillary incisor in two-jaw surgery. *J Oral Maxillofac Surg* 1984;42:402–404.
  44. Henry CH. Modified boley gauge for use as a reference plane in orthognathic surgery. *J Oral Maxillofac Surg* 1990;48:535–539.
  45. Schweska R, Engelke D, Zimmer B, Kubein-Meesenburg D. Positioning control of the upper incisors in orthognathic surgery: Preoperative planning with the Model Positioning Device and intraoperative application of the Sandwich Splint. *Eur J Orthod* 1991;13:367–371.
  46. Erickson KL, Bell WH, Goldsmith DH. Analytical model surgery. In: Bell WH (ed). *Modern Practice in Orthognathic and Reconstructive*

- tive Surgery. Philadelphia: Saunders, 1992: 154-205.
47. Schweska-Polly R, Roesse D, Kuhnt D, Hille K-H. Application of the model-positioning appliance for three-dimensional positioning of the maxilla in cast surgery. *Int J Adult Orthod Orthognath Surg* 1993;8:25-31.
  48. Turvey TA. Simultaneous mobilization of the maxilla and mandible: Surgical technique and results. *J Oral Maxillofac Surg* 1982; 40:96-99.
  49. Johnson DG. Intraoperative measurement of maxillary repositioning: An ancillary technique. *Oral Surg Oral Med Oral Pathol* 1985;60:266-268.
  50. Nishioka GJ, van Sickels JE. Modified external reference measurement technique for vertical positioning of the maxilla. *Oral Surg Oral Med Oral Pathol* 1987;64:22-23.
  51. Kahnberg K-E, Ridell A. Combined Le Fort I osteotomy and oblique sliding osteotomy of the mandibular rami. *J Craniomaxillofac Surg* 1988;16:151-156.
  52. Neubert J, Bitter K, Somsiri S. Refined intraoperative repositioning of the osteotomized maxilla in relation to the skull and TMJ. *J Craniomaxillofac Surg* 1988;16:8-12.
  53. Kahnberg K-E, Sunzel B, Astrand P. Planning and control of vertical dimension in Le Fort I osteotomies. *J Craniomaxillofac Surg* 1990; 18:267-270.
  54. Perkins SJ, Newhouse RF, Bach DE. A modified boley gauge for accurate measurement during maxillary osteotomies. *J Oral Maxillofac Surg* 1992;50:1018-1019.
  55. Van Sickels JE, Larsen AJ, Triplett RG. Predictability of maxillary surgery: A comparison of internal and external reference marks. *Oral Surg Oral Med Oral Pathol* 1986; 61:542-545.
  56. Stanchina R, Ellis III E, Gallo WJ, Fonseca RJ. A comparison of two measures for repositioning the maxilla during orthognathic surgery. *Int J Adult Orthod Orthognath Surg* 1988;3:149-154.
  57. Schweska R, Engelke D, Kubein-Meesenburg D, Luhr HG. Control of vertical position of the maxilla in orthognathic surgery: clinical application of the sandwich splint. *Int J Adult Orthod Orthognath Surg* 1990;5:133-136.
  58. Schweska R, Roesse D, Kuhnt D, Walloschek W. Splint for controlling vertical position in maxillary osteotomies. *J Clin Orthod* 1990; 24:427-431.
  59. Schweska R, Engelke D, Kubein-Meesenburg D. Condylar position control during maxillary surgery: the condylar positioning appliance and three-dimensional double splint method. *Int J Adult Orthod Orthognath Surg* 1990; 5:161-165.
  60. Schweska-Polly R, Engelke D, Kubein-Meesenburg D. Application of the condylar positioning appliance in mandibular sagittal split osteotomies with rigid skeletal fixation. *Int J Adult Orthod Orthognath Surg* 1992; 7:15-21.

