S-A-M®

MSF

Magnet  
Splitcastform

SAM®

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DF-Patent 3117506

MANUALE SISTEMA  
SPLIT-CAST MAGNETICO

*MANUAL MAGNET  
SPLITCAST FORMER SYSTEM*

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*RED: normal-sized former plate for primary base Ø 80 mm*

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*GREEN: extra large for a complete prosthesis Ø 100 mm*

## THE SAM® MAGNET-SPLITCAST-FORMER

The MSF system (magnet splitcast former system) is used to produce casts with a separating base and a magnet-secured primary/secondary base for

- Splitcast technology, to check correct assembly in the articulator.
- Comparison of centric impressions, for which the magnet can be removed from the base.
- Pinned casts with removable dies.
- Dental casts having soft tissue borders

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YELLOW: **small** base formers for smaller dental impressions or casts Ø 70 mm

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BLACK: **reciprocal** base plate for prior preparation of secondary cast base Ø 80 mm

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### PRODUCTION OF A SPLITCAST USING A DENTAL IMPRESSION

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*The yellow splitcast former is used in the following description. The other SAM® splitcast forms are to be used in the same way.*

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*The preparation of each form plate starts by placing the adhesion plate between the conical positioning pins in the centre of the form plate. The adhesion plate is secured by the magnet in this position. The rubber ring is then placed around the form plate.*

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*The cast and the splitcast former are poured with vacuum mixed dental gypsum.*

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*The filled magnet splitcast form is placed onto the upper poured impression. The base former should be distally lower than anterior.*

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
*Once the gypsum has set completely, the form plate is removed. The magnetic holding cup is placed onto the adhesion plate of the primary base. The gypsum surface of the primary base is saturated with SAM® UG-Sep (gypsum separator). After obvious saturation, pour of the excess and dry completely with compress air.*

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*To produce the secondary base, the rubber ring is filled with gypsum up to the upper edge. This is also done with the aid of a vibrator. Place a 10 CM square of small bubble wrap onto the poured gypsum to produce retention for articulator attachment.*

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*Once the gypsum has set completely, the rubber ring is removed. The split-cast can then be trimmed. The primary base should not be removed from the secondary base prior to trimming.*

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*Once the trimming has been completed, the bases are easily separated for the first time, rinsed with water, dried with compressed air and reassembled.*

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
### PRODUCTION OF A PIN CAST USING A SPLITCAST

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*The harmony between the length of the pin sleeves and the edge of the rubber ring is a unique feature of the SAM® AXIOPIN® system.*

*A small amount of vaseline placed on the tip of the AXIOPIN® after the pin sleeve has been set in place makes it easier to identify and make visible the pin tips on the under surface of the poured cast.*

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*Pinned dental casts are sprayed with SAM® UG-Sep, after saturation, dry completely with compressed air.*

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*The rubber ring is placed onto the form plate and a magnetic cup with a magnet or an adhesion plate is positioned, which is then filled to the upper rim with a type of gypsum with special flow properties (WhipMix „FlowStone FS” or Dentona „sockel-plaster GT 160”).*

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*To speed up cast production, the secondary base can be pre-made, using a reciprocal base former. This is done by using the black reciprocal splitcast former with convex ribs to obtain a secondary base, the ribs of which have been inserted into grooves, so that the pinned casts can be optimally positioned with many pins.*

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*The adhesion plate or a magnetic holding cup with a magnet are placed into the center of the magnetic base former and the higher, red rubber ring is used.*

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*The base form is filled with gypsum up to the edge.*

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*Bubble wrap is placed onto the gypsum base to ensure retention.*

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*Once the gypsum has hardened, only the reciprocal form plate is removed. The gypsum base is provided with a magnetic cup and insulated using gypsum separating agent (SAM® UG-Sep).*

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*On this prepared gypsum base is positioned e.g. a pin model...*

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*Unlike the „normal“ splitcast former (yellow, red, green), cast production using the reciprocal form plate (black) means producing a secondary base first, followed by the primary cast.*

*The gypsum can be poured directly to the prepared, insulated gypsum split-cast plate.*

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*The boxing aid is used for the production of casts having border moulded soft tissues and a magnet-secured primary and secondary base.*

1 Basi split-cast XL  
*two splitcast formers*

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*The impression is placed onto the base plate, using three supports made of moldable elastic modelling clay. The front tooth area of the upper impression is lowered, while for the lower impression, the molar area is the lowest part.*

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*The black sleeve is placed onto the base plate. The distance gauge is placed onto the sleeve.*

*Optimal cast thickness is guaranteed when the pin touches the highest arch of the impression; for the upper impression, this is the apex of the palate, while for the lower impression it is the edge on the lingual side.*

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*The black plastic sleeve is sealed flush with the green clamping piece (arrow) and thus fixed on the base plate.*

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*Thinly mixed alginate is now used to fill up the space at least up to the functional margin.*

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*Once the alginate has  
hardened, the black  
sleeve is removed and the  
verfilled alginate is  
trimmed up to the  
required extension of the  
functional margins.*

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*The trimming of the alginate determines the extension of border moulded soft tissues*

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*The sleeve is then placed back around the alginate bed, lowered onto the base plate and sealed with the clamping piece.*

*The idea is to obtain functional margins for the overall prosthetic cast in which the nested functional impression border defines the limits without further processing.*

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*The second splitcast form plate is equipped with an adhesion plate, layered with gypsum and inserted into the plastic sleeve. The splitcast form is lowered into the sleeve until it is flush with the edge.*

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*Any excess gypsum is immediately removed.*

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*Once the gypsum has hardened, the magnetic splitcast form is removed and any remaining gypsum is removed from the adhesion plate. The gypsum surface of the primary base is saturated with gypsum separating agent (SAM® UG-Sep) and dried completely with compressed air.*

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*The magnet is placed onto the adhesion plate, using the magnetic holding cup.*

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*The secondary base is cast by using a vibrator to fill the cast with gypsum up to the rim of the sleeve.*

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*Bubble wrap is placed onto the liquid gypsum to ensure retention.*

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*Once the gypsum has hardened, the base plate, sleeve and alginate bed are taken apart.*

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*First the base edges are trimmed, followed by the underside, after which the base plates of the cast are separated from each other **for the first time**, rinsed with water and dried with compressed air.*

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*Once the upper and lower jaw casts have been assembled in the articulator, the splitcast may be tested.*

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
*To do this, the incisal pin is lowered and the magnet lifted out of the base, using a magnet remover.*

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*The secondary base in the upper section of the articulator is lowered onto the upper jaw model with the centric lock closed.*

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*The cast has been correctly assembled when there is a good fit between the primary and secondary bases when the articulator is closed.*

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*Where there is a gap between the primary and the secondary base, this indicates an incorrect assembly. This is also true for the comparison of different records.*

AXIOGRAPH® • AXIOSPLIT® • AXIOQUICK®  
AXIOSAW® • AXIOPIN® • AXIODRILL® • AXIOSIM®  
AXIOCOMP® • AXIOWAX® • AXIOTRON® • SAM®

sono marchi registrati di SAM® Präzisionstechnik GmbH  
*are registered trade marks of the SAM® Präzisionstechnik GmbH*

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