SAM

Magnet
Splitcastform

SAM[®]

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DE-Patent 3117 506

**MAGNET SPLITCAST
FORMER**

INHALT

CONTENTS

	Magnet-Splitcast-Former	3
A -	Herstellung eines Splitcast-Modelles mit Abformlöffel <i>Production of a splitcast using a dental impression tray</i>	10
B -	Herstellung eines Pin-Modelles mit Splitcast <i>Production of a pin cast using a splitcast</i>	18
C -	Sockelformer, reziprok <i>Reciprocal base former</i>	22
D -	Modelle zur Erhaltung der Funktionsränder <i>Models to preserve the functional edges</i>	30
E -	Die Splitcastkontrolle <i>Splitcast testing</i>	46
	Pflege Tipps <i>Maintenance tips</i>	51

MAGNET-SPLITCAST-FORMER

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MAGNET-SPLITCAST-FORMER

The MSF system (magnet splitcast former system) is used to produce casts with a separating base and a magnetsecured primary/secondary base:

- The splitcast technique is used to check the model assembly.*
- By removing the magnet, centric registers can be compared.*
- Production of pin models with splitcast base.*
- Making a Precise Function Wheel Model with Splitcast.*

Alle Splitcast Former Kits enthalten:

- 6 Magnete
- 6 Haftplatten
- 6 Magnettöpfe
- 1 Magnetheber

All splitcast former kits include:

- 6 Magnets*
- 6 Adhesion plates*
- 6 Magnet holding cups*
- 1 Magnet remover*

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YELLOW

Size S - Ø 70 mm

Small base formers for offices without trimmer or for small jaws

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BLACK

Size M - Ø 80 mm

Reciprocal former plate for prior preparation of secondary cast base

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GREEN

Size XL - Ø 100 mm

Former plate for total prosthetics and large jaws

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PRODUCTION OF A SPLITCAST-MODEL WITH IMPRESSION TRAY

Magnet Splitcast Former Produkte
Magnet splitcast former products

Splitcastformer
Splitcast former

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The preparation of each form plate starts by placing the adhesion plate between the conical positioning pins in the centre of the form plate. The adhesion plate is secured by the magnet. In this position the rubber ring is then placed around the form plate.

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The cast and the splitcast former are poured with vacuum mixed dental gypsum.

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The poured impression is centered on the filled magnet-splitcast-former set. It is recommended that the model be turned over when the plaster appears dull of hardening. The base former should sit distally deeper than anteriorly.

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After the gypsum has set completely, the form plate can be removed. The magnetic holding cup is placed onto the adhesion plate of the primary base. The gypsum surface of the primary base is saturated with SAM® UG-Sep (gypsum separator). After obvious saturation, pour out the excess and dry completely with compress air.

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To create the secondary base, the rubber ring is filled with plaster up to the upper edge. To form retentions, a piece of bubble wrap is pressed into the wet plaster.

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Once the gypsum has set completely, the rubber ring have to be removed and the splitcast can then be trimmed. The primary base should not be removed from the secondary base prior to trimming.

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Only after trimming may the base be separated, rinsed with water, dried with compressed air and reassembled.

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
PRODUCTION OF A PIN CAST WITH SPLITCAST FORM

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A special feature of the SAM® AXIOPIN® system is the length of the pins and sleeves. These are matched to the edge height of the rubber ring.

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The manual and more detailed information to AXIOPIN® Pins and Sleeves and the pin drill AXIODRILL® are available from SAM or as PDF for download:

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Pinned dental casts are sprayed with SAM® UG-Sep. After saturation, dry it completely with compressed air.

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The rubber ring and an adhesive plate are placed on the mould plate. Then the ring is filled with base plaster up to the upper edge. Otherwise, plaster can be applied with a brush or probe at the gear rim between the sleeves. The gear rim is then placed in the base former.

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To speed up model fabrication, a secondary base can be pre-fabricated with excess plaster. SAM offers the black reciprocal splitcast former with convex ribs for this purpose. This secondary socket, which has the ribs as indentations, simplifies the positioning of a dental arch with many pins.

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The adhesive plate is inserted in the center of the base former. The higher, red rubber ring is used for the black forming plate.

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The base former is filled up to the edge with plaster. It is best to vibrate the plaster in so that the groove retention of the adhesive plate is completely filled.

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For retention formation, the nubs of the bubble wrap are pressed into the wet plaster base..

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Once the gypsum has hardened, only the reciprocal form plate is removed. The remaining plaster base is insulated with plaster release agent (SAM® UG-Sep) and fitted with a magnetic pot.

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For example, a plaster plate with many pins can be optimally positioned on the resulting plasterboard.

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Note: The molded sockets of red and black base formers do not match.

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In contrast to the yellow and red socket former, the reciprocal socket former (black and green) is used to create a secondary socket first and then the primary model. The plaster can be applied directly to the prepared, isolated plaster splitcast plate.

Tip: The reciprocal base former is ideal for recycling excess plaster.

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FUNCTIONAL MODEL WITH SPLITCAST BASE

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The single box aid is used to produce models with a functional rim and a splitcast base.

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The functional impression is placed on the base plate with a three-point support made of permanently elastic plasticine rubber.

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The black collar is placed on the base plate. The distance gauge is placed on the collar. The optimum model thickness is achieved when the pin touches the highest curvature of the impression. For the upper jaw, this is the palatal margin, for the lower jaw it is the lingual margin.

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The black plastic sleeve is going to be sealed precisely using the green clamping piece (arrow). Apply it from the top, then pull it down at the bottom plate.

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The space is filled up to the functional edge with thinly mixed alginate. Please observe the setting time specified by the manufacturer.

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After the alginate has hardened, remove the black sleeve, than trim up the overfilled alginate up to the required extension of the functional margin.

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The plastic collar is then placed around the alginate bed again, lowered to the base plate, closed with the clamping sleeve and filled with super hard plaster.

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The second splitcast form plate is equipped with an adhesive plate. The retentions of the adhesive plate should be covered with plaster to avoid bubble formation. The plate is then inserted into the plastic sleeve and lowered to the edge of the plate under vibration.

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The excess plaster should be removed immediately.

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After the plaster has hardened, only the magnet splitcast form is lifted off and the adhesive plate is freed from plaster residues. The plaster surface of the primary base is isolated with the plaster release agent SAM® UG-Sep and blown dry after 30 seconds.

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The magnet with magnet cup is placed on the holding plate.

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During production of the secondary base, the plaster is vibrated up to the edge of the plastic collar.

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For retention formation, bubble wrap is pressed into the wet plaster.

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Once the gypsum has hardened, the base plate, sleeve and alginate bed are taken apart.

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First the base edges and then the underside are trimmed. Then the base plates of the model are separated for the first time, rinsed with water and dried with compressed air.

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For checking the model assembly with centric registration.

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Once the upper and lower jaw casts have been assembled in the articulator, the splitcast may be tested.

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
The incisal pin is lowered and the magnet removed from the base with a magnetic remover.

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The secondary base in the upper section of the articulator is lowered onto the upper jaw model with the centric lock closed.

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The cast has been correctly assembled if there is a good fit between the primary and secondary bases when the articulator is closed.

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Where there is a gap between the primary and the secondary base, this indicates an incorrect assembly. This is also true for the comparison of different records.

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MAINTENANCE TIPS

The splitcast former plates and rubber rings can be easily cleaned with a toothbrush and running water. For persistent gypsum encrustation this can be solved beforehand by soaking in water. Vaseline can be used to prevent the rubber rings from erosion in the long run.

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*Universal release agent
 against plaster.
 To isolate the plaster surface,
 spray it with UG-Sep and wait
 30 seconds, then dry the plaster
 with compressed air.*

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*Adhesion plates with grooved
retention
(50 pieces)
for magnetic
splitcast former*

REF MSF 120

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*Magnet holding cups
(50 pieces)
for MSF magnets*

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*Magnets for MSF
(20 pieces)
with MSF magnet remover*

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Magnet splitcast former combination kit

*Includes all four splitcast former
yellow S (Ø 70 mm), red M (Ø 80 mm), black reciprocal M (Ø 80 mm)
and green XL (Ø 110 mm), 4 rubber rings, 1 magnet remover,
6 magnets and 6 magnet cups, 6 adhesion plates and 1 UG-sep spray.*

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Of course, the splitcast former can also be ordered separately (see page 6-9).

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The SAM® AXIOPIN® pin system is used for the production of pinned models for the removal of individual tooth model dies.

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parallel dowel pin system, with metal pins and yellow plastic sleeves, 1,6 mm drill hole size, dowel pins are attached to the dental cast with cyanoacrylate glue

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Horizontal pin drilling machine with a view of the drill tip and optimum positioning of the drill hole according to the saw cut drawings, with drilling depth gauge and hexagon key.

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are registered trade marks of the SAM® Präzisionstechnik GmbH
sind eingetragene Warenzeichen der SAM® Präzisionstechnik GmbH
SAM® Präzisionstechnik GmbH