# SAM 2P

HANDBUCH  
*MANUAL*

# Inhalt

## Contents

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
*The **Upper Member** of the articulator consists of **Condylar Housings** posteriorly and Upper **Incisal Guidance Table** anteriorly. Together, they simulate the TMJ and upper anterior teeth guidance.*

*The moving **Lower member** of the articulator with **Condylar Elements** posteriorly and the **Incisal Pin** anteriorly, together, they simulate the mandibular condylar movement and lower anterior teeth guidance.*

*This, therefore, provides for a very logical and natural laboratory simulation of the patient movements which are very easy to understand and reproduce.*

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*The red dot in the center of the lower member condylar element corresponds to a hinging axis rotation center of the mandibular condyles. Also known as "Hinge Axis".*

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*Above the red dotted/dashed line is the axis-orbital plane (AOP)*

*The upper anterior incisal table surface along with the posterior condylar element centers of rotation are on the same horizontal plane.*

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
*The lower incisal pin in contact with the upper incisal guidance table is used to simulate the relationship of upper and lower incisal edges of the patient's teeth. Thus promoting normal functional patterns.*

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

*All current articulators and instruments made by SAM are based on the height „P“ (prosthetic height = 121 mm between upper and lower members). The intercondylar distance is 110 mm. A particular feature of the system is the SAM<sup>®</sup> 15 mm distal offset arrangement of the mounting plates in the articulator. This creates a more normal patient simulation of mounted models in the SAM articulator.*

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
*The SAM<sup>®</sup> 2P and all articulators made by SAM are ARCON Instruments. To simulate the patient's authentic TMJ movements, the lower articulator member condylar ball moves in the upper member condylar housing and this movement corresponds to the TMJ anatomical model. This natural anatomical movement is logically simulated and understandable.*

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*The SAM 2P is a patient simulator which allows the end-user to reproduce anatomic patient mandibular movements for examination, diagnosis, and therapy.*

Aus dem digitalen Artikulator AXIOCOMP  
*Excerpt from the digital articulator AXIOCOMP*

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
*The articulator axis orbital (Frankfort Horizontal) plane and all horizontal planes of the articulator are parallel with the table top. Centric relation mounted casts are anatomically related to the patient's cranium with enough space for all types of casts, including those with split casts.*

*The SAM 2P articulator has the same functional dimensions as all other SAM articulators (Arcon principle).*

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*Single articulators and small kits are delivered in a translucent carrying case, larger kits are in the yellow carrying case.*

*In addition, articulator mounted casts can be transported in the carrying cases along with the articulators.*

*It is recommended that all SAM instrumentation, parts, etc. should be placed in the designated areas in the foam to ensure safe transport or dust-free storage.*

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*Thumb Screw for the Inclination of the Condylar Housing*

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
Haftplatte  
*Adhesion Plate*


Magnet-Montageplatte  
*Magnetic Mounting Plate*

# Artikulator Unterteil *Articulator Lower Member*

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
Schraubmontageplatte  
*Screwtype Mounting Plate*

Magnet mit Magnettof  
*Magnet with MPS Magnetic Holding Cup*

AXIOSPLIT® System  
*AXIOSPLIT® System*

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*Spring Loaded Screw of Centric Lock*

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
Draufsicht / *Top View*

*Centric is unlocked, if spring loaded screw is screwed outwards.*

*Turn the Screw of the Centric Lock **Inwards** to Lock the Centric.*

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Draufsicht / *Top View*

*Centric is Locked, if Spring Loaded Screw is in touch with Condylar Ball.*

*Turn out the Screw of the Centric Lock **Outwards** to Unlock the Centric.*

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
## Der erste Digitale SAM Artikulator *First Digital SAM Articulator*

vom Erfinder / by inventor of

AXIOGRAPH®, AXIOTRON® - AXIOQUICK® Recorder

**Unentbehrlich in Didaktik und Lehre!!**  
*Indispensable in Didactics and Teaching*

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*Lower jaw movements are the result of rotation and translation in the jaw joints. In the SAM articulator, the horizontal condylar (articular path) inclination (HCI) can be adjusted, for median movement a skew angle (Bennett) can be set. The individual setting according to measured data is supplemented by inserts with curvature (HCI & Bennett) in the modular SAM system.*

*With the SAM transfer system, the anatomic correct model transfer in the articulator is ensured, a prerequisite for the rotation and translation of the jaw joints.*

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#### *Green (included)*

*For Normal initially curved Sideshift Movements*

0,5 mm:	54°
3,0 mm:	15°
10,0 mm:	7,5°

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
#### *Red (optional)*

*For Immediate Sideshift Movements*

0,5 mm:	78°
3,0 mm:	38°
10,0 mm:	13°

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
*The Bennet Angle can be changed from 0°- 25° when using the straight Bennett Guides. If color-coded Bennett Guides are used, the Bennet angle can be varied from 0 - 10 °.*

# Einstellung des Bennettwinkels

## Setting the Bennett Guidance Angle

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
*The figure beside shows the effect of a curved Bennett Insert Guide with an initial inclination angle at 0.5 mm of the mediotrusion movement and 0° angle adjustment. Here, the first one to two millimeters are dominant.*

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
*Bennett Angle  
Bennett Angle scale and Adjustment Scale (scale in degrees); additional to preset Bennett Angle of the Green Insert*

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## *Changing Bennett Guidance Inserts*

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*Open the Knurled Screw of the Bennett Angle Guide and unscrew it completely.*

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*The Bennett Angle Guide Insert can be removed by light pressure on the threaded pin.*

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*Pay attention to the lateral correctness of the Bennett Angle Guide Inserts and insert it into the provided hole. (Arrow).*

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*Tighten the Knurled Screw of the Bennett Angle Guide Inserts and adjust the Bennett Angle (i. e. 5°) before final closing.*

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*Average Condylar Pathway Inclination is 35°. The average value relates skull to the Frankfort Horizontal Plane (Axis Orbital Plane). With Complete Dentures an average value of 20° is recommended.*

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*Setting the Condylar Pathway Inclination. The adjustment of the inclination is performed by loosening the posterior Thumb screw on the rear of the upper member of the articulator.*

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*The inclination of the Condylar Pathway is adjustable from 15° to 65°.*

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*For the simulation of the Condylar Pathway in the SAM 2P there are three different condylar housings with different Condylar Curvatures. These are numbered and interchangeable.*

*They differ in the curvature of the Protrusion pathways from each other. The convex inner guidance of the Condylar Housing is tracked by the Condylar Ball. This simulates the Protrusive and Lateral Movements of the lower jaw of the patient in the articulator.*

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



*To change the Condylar Housing and thus the curvature, the Knurled Screw for the Condylar Pathway locking is first loosened. The articulator upper member is separated from the lower member.*

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*The Knurled Locking Screw will be removed completely along with washer and the brass pressure piece.*

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*The Locking Screws will be unscrewed with a flat head screwdriver.*

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
*The Condylar Housing can now be removed by pushing and pulling. Please check for the numbered pairs which indicate the corresponding curvature.*

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*Rotate the Condylar Housing as it is being inserted into the opening on the correct side of the upper member as it snaps into place.*

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*If the Condylar Housing cannot be inserted, the bolt (arrow), which is secured with a spring, must be pressed downwards with an allen wrench. This allows the housing to be inserted without any problems. (See p. 38)*

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*The brass pressing piece is placed with the beveled side of the condylar axis. After this, insert the washer and insert the thumb screw. Holding this in place with the locking screw.*

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
*Finally the position of the housing must be secured again with the condylar housing locking screw*

*We recommend an average setting of 35° for the Condylar Pathway Inclination.*

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*In order to adjust comfortably selected positions during protrusive and lateral motions in the articulator, color coded protrusions inserts are used for the SAM 2P.*

*These are available in millimeter increments from 1 to 6.*

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*Set the White Straight Bennett Angle to 0°, since the protrusion movement corresponds to a straight forward movement. For a protrusion or retrusion movement the straight Bennett inserts must be used and the Centric Locking screw must be open.*

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*When the protrusion clip is inserted, the ridged surface must show outwards. The smooth side represents the pathway of protrusion, here it is green = 3 mm.*

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*The insert is hooked up onto the inner wall of the housing and pressed down with light pressure until the insert snaps.*

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
*Using the protrusion insert, the Condylar Ball can be moved forward and simulates a lower jaw movement straight forward.*

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
*Using the same protrusion clips on both sides, a protrusion can be adjusted up to 6 mm.*

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*If the inserts are used only on one side or different clips are inserted, the position must be defined in the median movement.*

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*The upper member feet prevent the rotation and must be removed with a screwdriver.*

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*The Centric Lock on both sides must be completely Locked.*

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*The Locking of the Condylar housings must be loosened on both sites completely.*

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*Push the cup of the vertical support rod onto the thumb screw of the incisal table.*

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*The Tilt Support Rod supports the upper member of the articulator in the open position. It can be opened by 180°.*

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*This improves the view of the occlusion, and makes it easier to work on the model (e.g wax up).*

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*Occasionally apply silicone spray to maintain the Tilt Support Rod and the inner O ring (ART 525002).*

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*A special feature of the SAM® Plate System is the 15 mm offset arrangement of the mounting plates in the articulator which takes into account the normal lower jaw distal relationship to the upper jaw. By this distal offset of mounting plates position, the mounting of the lower jaw model is easier.*

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*The different SAM Mounting Plate Systems are available for all SAM articulators and equipment. All SAM instruments can be converted without difficulty.*

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#### Mounting Plate with Screw Type System

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*The yellow mounting plate is secured with a Thumb Screw and this provides maximum retention for assembly of models with small amounts of plaster. The mounting plates are intended for single use and for storage and easy reattachment of the mounted models.*

*The low cost of the mounting plates indicates this is the procedure of choice. Not to mention the saving of the costly time to try and recycle the plates. The positioning of dental cast in the articulator is carried out very precisely by means of two precision alignment pins.*

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#### Benefit:

*the patient model remains from beginning to end in the same articulator*

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*The MPS Magnetic Mounting Plate System is an alternative to the historical screw mounting plate. This Magnetic Plate System allows for the quick removal and replacement of the models in the articulator. The conversion of screw type to magnetic system can be done easily and quickly with the appropriate kits (MSP 101K and 101K MPS).*

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*The AXIOSPLIT® system is used for trouble-free and accurate exchange of models between different articulators (for example, old and new devices), in practice and laboratory. So the models need only to be sent. The articulators stay in place.*

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*A Master Check ensures the synchronizing with all SAM devices. The mounted models can be interchanged without loss of fitting accuracy in all SAM® devices. You will initially need the AXIOSPLIT Starter Kit with Master Check, norm incisal pin and related matrix holder and an articulator kit for each additional instrument. The exact procedure and more information is described in our AXIOSPLIT-Manual.*

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**Benefit:** Highest precision and your own synchronizing

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## Screw Type Mounting Plates Attachment

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*The Screw Type Mounting Plate has a round hole behind the threaded nut and an elongated slot anteriorly to assure precise plate positioning on all SAM instrumentation set to the precision mounting plate pins. The elongated hole shows to the front.*

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
*With one hand the plate is held in position and with the second hand thumb and fore finger the thumb screw is tightened to contact. It is further tightened firmly only 1/4 of a turn. tightened with the knurled screw on the other side.*

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*The mounting plate is fully seated when it sets firmly on the border areas without gap.*

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*Take note: The mounting plate is screwed firmly to the articulator, the mounted model may not be set easily and quickly removed on and off as like with the MPS system.*

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*Note: The mounting plate is disposable and should not be used more than once since the difficulty of removing set gypsum and the cleaning of the mounting plate normally will destroy plate and destroy it's accuracy.*

# MPS-Magnet Montageplatten-Befestigung

## MPS-Magnetic Mounting Plates Attachment

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*The adhesion plate is set with the rounded edge downwards into the provided hole.*

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
*The MPS-mounting plate is attached with the round side containing the slot anteriorly onto the positioning pins.*

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
*Take note: The slightly rounded and shiny side shows downwards when it is inserted into the mounting plate (the rounded edge side shows space for plaster).*

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*To reproduce or build tooth guidances there should to be an incisal table. It must be able to produce the inclination and direction of the tooth guides.*

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
*To simulate the anterior teeth and canine guidance authentically, the incisal table is always attached to the upper part of the articulator.*

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*The incisal table is located as a guiding surface on the top member of the articulator. The flat plane incisal table is used to secure the vertical dimension and to support the incisal pin and can be individualized with plastic via the lateral wings.*

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*The incisal table is not pushed completely into the articulator top until it stops until the anterior moveable section has been completely raised.*

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*Set incisal table on steepest inclination 80° and the rotation to 0°*

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*The incisal table is seated fully in the slot on the upper member of the articulator until the tip of the incisal pin simultaneously touches the inclined and the horizontal surface of the incisal table. Then the knurled screw is tightened.*

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*In order to check the correct position of the incisal table, the incisal table can be rotated anywhere from 0° to 60° around its own axis when the articulator is closed.*

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*The incisal pin should be located in the lower part of the SAM articulator.*

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
*The zero mark on the scale in the middle of the incisal pin is normally the starting position for the dental laboratory.*

*The incisal pin can be fixed in the restorations bite height. The incisal pin can be raised to 11 mm or lowered up to 6 mm.*

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*Due to the articulator construction the desired value must be set twice at the incisal pin while raising the occlusion, so e.g. in an increase of 2.5 mm in the occlusion the adjustment at the incisal pin is +5 mm.*

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*When lowering of the bite to 1 mm, the incisal pin should be lowered 2 mm.*

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*The micro adjustable incisal pin can be adjusted accurately and within 0.05 mm range. When lateral knurled screw is opened the pin can be finely adjusted on the upper black scale of incisal continuously.*

*This can be finely adjusted in the micrometer range.*

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*To **Raise**, turn **Clockwise**.  
Rotate **Counterclockwise** to **Lower**.*

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*Open the lateral knurled screw slightly (about 1/8 turn).*

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*Rotate the black micro adjustment screw clockwise until the horizontal mark on the scale is in the center of the pin and shows "2" with the vertical mark on the upper black scale "3".*

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*The position is fixed with the lateral knurled screw.*

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

*The axis orbital plane in the patient corresponds in the articulator to the line between the bearing surface of the incisal table to the incisal pin (1) and the centers of the condylar balls (hinge axis) (2).*

*Cranium-related mounted models in the articulator are an essential part of diagnosis and therapy. The correct position of the maxillary model to the reference plane is an absolute prerequisite for model mounting with a centric registration, but also for setting condylar pathways and even when using mean values.*

*The Frankfort Horizontal (FH) Plane represents a reference plane for the patient and the articulator. The reference points of the Frankfort horizontal (Porion and orbitale) are identical to the axis-orbital plane (AOE).*


*The scissors mechanism of the SAM<sup>®</sup> anatomic facebow creates and allow for an optimal positioning of the transfer fork (from a functional and aesthetic point of view).*

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*Please find the exact procedure and the correct technique to attach the facebow and the proper maintenance and care in our AXIOQUICK facebow manual. This manual can be provided on request by SAM or can be downloaded as PDF file.*

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*The transfer stand AX is used for an easy and fast maxilla mounting in all SAM articulators and laboratory cast mounting instruments. After the transfer the anatomic facebow can remain in the treatment room, only the transfer fork assembly is needed for the model transfer. All versions of the SAM transfer fork assemblies are compatible with the transfer stand AX. With the transfer stand the upper jaw position of the patient can be quickly and precisely transferred to the articulator.*

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*A more detailed manual for model mounting can be provided on request by SAM or can be downloaded as PDF file.*

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*The occlusal plane indicator set (ART 385) with a plane and a curved plate with Curve of Spee set can be used for the average-maxillary model mounting or for a prosthetic set up.*

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*The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm.*

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
*The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm The angle scale can be adjusted between 0° - 45°, we recommend that the inclination is usually between 0° - 15°.*

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*The template is used in the production of complete dentures and as an aid for a setting up of the maxillary teeth. It provides a Curve of Spee. The marks on the template indicate auxiliary lines as the middle of jaw ridge. In the loosened state, the cap may also adjust the inclination angle of the remaining teeth can be determined as occlusion plane.*

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## *Cleaning and Maintenance of the Articulator*

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*Clean the centric lock and the condylar housing with cleaning solvent to eliminate residues of grease, metal and dust.*

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*The cleaning solvent and the dirt is wiped off with a dry cloth.*

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*Lubricate the moving parts with silicone grease (ART 193)*

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*Then spread the lubricant while moving in the centric lock.*

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*Every now and then, the axis of the condylar housings should be freed from dirt and grease residues, so that the housing remains smooth-running.*

*Remove successively the thumb screws (A) for the inclination of the condylar housings, the washer (B) and the brass pressure piece (C).*

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*The locking screws can be unscrewed with a screwdriver.*

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*The condylar housing can now be removed by pushing and pulling.*

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*All metal parts are cleaned with cleaning solvent and freed from debris.*

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*After cleaning, all moving metal parts should be lubricated with a silicone grease e.g. like this aluminium shaft.*

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*While rotating the condylar housing it is inserted into the shaft hole, and snaps into place.*

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*If the condylar housing does not slip in, press the securing bolt (secured with a spring) down for example with an allen wrench and slide the housing in.*

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*All repairs, maintenance, adjustments and engravings or laser inscriptions on your device can be carried out professionally and quickly by our service team.*

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*For cleaning and maintaining your articulator we recommend our field checking kit (ART 187)*

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*If you need additional mounting plates or if you have lost a screw, a nut, a washer or any other item from your SAM equipment, please get in contact with us or your local dealer. They have a small stock of most common parts. Don't hesitate to contact them!*

*You can also download our catalogues for products, spare parts and consumables from our download-page:*

[http://www.sam-dental.de/pages/de\\_Catalogues.html](http://www.sam-dental.de/pages/de_Catalogues.html)

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SAM® Präzisionstechnik GmbH ist
*SAM® Präzisionstechnik GmbH is
**SAM Präzisionstechnik GmbH**