SAM 3

HANDBUCH *MANUAL*

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
*The **Upper Member** of the articulator consists of **Condylar Housings** posteriorly and Upper **Incisal Guidance Table** anteriorly. Together, they simulate the TMJ and upper anterior teeth guidance.*

*The moving **Lower member** of the articulator with **Condylar Elements** posteriorly and the **Incisal Pin** anteriorly, together, they simulate the mandibular condylar movement and lower anterior teeth guidance.*

This, therefore, provides for a very logical and natural laboratory simulation of the patient movements which are very easy to understand and reproduce.

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The red dot in the center of the lower member condylar element corresponds to a hinging axis rotation center of the mandibular condyles. Also known as "Hinge Axis".

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
Above the red dotted/dashed line is the axis-orbital plane (AOP)

The upper anterior incisal table surface along with the posterior condylar element centers of rotation are on the same horizontal plane.

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
The lower incisal pin in contact with the upper incisal guidance table is used to simulate the relationship of upper and lower incisal edges of the patient's teeth. Thus promoting normal functional patterns.

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
All current articulators and instruments made by SAM are based on the height „P“ (prosthetic height = 121 mm between upper and lower members). The intercondylar distance is 110 mm. A particular feature of the system is the SAM[®] 15 mm distal offset arrangement of the mounting plates in the articulator. This creates a more normal patient simulation of mounted models in the SAM articulator.

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
The SAM[®] 3 and all articulators made by SAM are ARCON Instruments. To simulate the patient's authentic TMJ movements, the lower articulator member condylar ball moves in the upper member condylar housing and this movement corresponds to the TMJ anatomical model. This natural anatomical movement is logically simulated and understandable.

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The SAM 3 is a patient simulator which allows the end-user to reproduce anatomic patient mandibular movements for examination, diagnosis, and therapy.

Aus dem digitalen Artikulator AXIOCOMP
Excerpt from the digital articulator AXIOCOMP

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
The articulator axis orbital (Frankfort Horizontal) plane and all horizontal planes of the articulator are parallel with the table top. Centric relation mounted casts are anatomically related to the patient's cranium with enough space for all types of casts, including those with split casts.

The SAM 3 articulator has the same functional dimensions as all other SAM articulators (Arcon principle).

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
Single articulators and small kits are delivered in a translucent carrying case, larger kits are in the yellow carrying case.

In addition, articulator mounted casts can be transported in the carrying cases along with the articulators.

It is recommended that all SAM instrumentation, parts, etc. should be placed in the designated areas in the foam to ensure safe transport or dust-free storage.

Articulator Upper Member


Optional: Protrusions-
schrauben statt Amboss

*Optional: Adjustable
Protrusion Screws*

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
*Thumb Screw for the
Inclination of the
Condylar Housing*

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
Optional einstellbarer Inzisaltisch
Optional Adjustable Incisal Table


Articulator Lower Member


Transferbogen-Ohrolive-Montagestift
Face Bow Earpiece Mounting Pin

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


AXIOSPLIT® System
AXIOSPLIT® System

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

Draufsicht / *Top View*

Centric is Unlocked, if Cone Ring is in touch with head of post.

*Push Cone Ring of the Centric Lock **Inwards** to Lock the Centric.*

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
Draufsicht / *Top View*

Centric is Locked, if conus ring is in touch with Condylar Housing.

*Push the Cone Ring of the Centric Lock **Outwards** to Unlock the Centric.*

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Condylar Housing Covers are normally removed during:

- Wax-up-Technik
Wax-up Technique
- Wechsel der Bennett-Führungen
Changing of the Bennett Guidance Inserts
- Wechsel der Kondylarbahn-Einsätze
Changing of the Condylar Pathway Inserts
- Reinigung des Kondylengehäuses
Cleaning the Condylar Housing

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The Condylar Housing Cover is used to keep Upper and Lower Members of the articulator in contact while Centric Lock is OPEN.

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To remove the Condylar Housing Covers, the Thumb Screw needs to be OPEN completely.

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Note: The left and right side covers are different.

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
Lower jaw movements are the result of rotation and translation in the jaw joints. In the SAM articulator, the horizontal condylar (articular path) inclination (HCI) can be adjusted, for median movement a skew angle (Bennett) can be set. The individual setting according to measured data is supplemented by inserts with curvature (HCI & Bennett) in the modular SAM system.


With the SAM transfer system, the anatomic correct model transfer in the articulator is ensured, a prerequisite for the rotation and translation of the jaw joints.

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Green metal (included)

For Normal initially curved Sideshift Movements

0,5 mm:	54°
3,0 mm:	15°
10,0 mm:	7,5°

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Red (optional)

For Immediate Sideshift Movements

0,5 mm:	78°
3,0 mm:	38°
10,0 mm:	13°

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
The bennet angle can be changed from 0°- 25° when using the straight Bennett guides .If color-coded Bennett guides are used, the Bennett angle can be varied from 0 - 10 °.

Setting the bennett guidance angle

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
The figure beside shows the effect of a curved Bennett Insert Guide with an initial inclination angle at 0.5 mm of the mediotrusion movement and 0° angle adjustment. Here, the first one to two millimeters are dominant.

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
*Bennett Angle
Bennett Angle Scale and Adjustment Tab (scale in degrees); additional to Breset Bennett Angle of the Green Insert*

Changing Bennett Guidance Inserts

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
Open the Knurled Screw of the Bennett Angle Guide and unscrew it completely.

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The Bennett Angle Guide Insert can be removed by light pressure on the threaded pin.

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Pay attention to the lateral correctness of the Bennett Angle Guide Inserts and insert it into the provided hole. (Arrow).

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
Tighten the Knurled Screw of the Bennett Angle Guide Inserts and adjust the Bennet Angle (i. e. 5°) before final closing.

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
Average Condylar Pathway Inclination is 35°. The average value relates skull to the Frankfort Horizontal Plane (Axis Orbital Plane). With Complete Dentures an average value of 20° is recommended.

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Setting the Condylar Pathway Inclination. The adjustment of the inclination is performed by loosening the posterior Thumb screw on the rear of the upper member of the articulator.

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The inclination of the Condylar Pathway is adjustable from 15° to 65°.

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
Three different Condylar Pathway Inserts are available for simulating the condylar pathway in the SAM 3. They are color-coded and interchangeable. They differ in the curvature of the protrusion pathway from each other.


The convex inner guidance of the curvatures in the Condylar Housing Inserts are tracked by the Condylar Ball. This simulates the Protrusive and Lateral Movements of the lower jaw of the patient in the articulator.

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

- Curvature 1:** HCI* 35°
- Curvature 2:** HCI* 30°
- Curvature 3:** HCI* 30°

* Horizontal Condylar Inclination

Changing the Curvature Pathway Inserts

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*The Bennett Angle is set to 0°.
Separate articulator Upper member from the Lower Member.*

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The Curvature Insert Locking Screw on the front side of the Condylar Housing is loosened by rotating it counterclockwise.

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
The curvature insert can now be removed. Select the correct right and left curvature pair. The correct inserts can be identified by markers (R1 + / L1 +) just under the screw on the inside of the insert. The following applies: the square cut in the curvature next to the metal thread is positioned towards the Bennett Guidance Insert.

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The position of the curvature insert is checked by slightly pressing. Once seated, the locking screw is inserted and tightened clockwise. The process is repeated on the opposite side. The articulator Lower Member can then be attached the Upper Member

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
The SAM 3 can be retrofitted with adjustable Protrusion Screws.

The Bennet Angle is set to 0°, since the protrusion movement corresponds to a straight forward movement.

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For a Protrusion or Retrusion movement the straight Bennett inserts must be used and the Centric Locking must be open.

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Protrusion can be adjusted up to 6 mm. Turn the SMALL Knurled Protrusion Screw counterclockwise. One line equals 1 mm. One rotation corresponds to 0.8 mm.

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
The initial position 0 is reached when the Centric Locking can be closed. The starting position must be taken before switching to the retrusion setting and vice versa.

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

Retrusion: The large Knurled of the Protrusion Screw is rotated counterclockwise during the retrusion. The lower jaw can be retruded up to 2 mm.

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*The Protrusion can be adjusted up to 6 mm by turning the small knurled screws counterclockwise on both sides.
If the screw is operated only on one side, the position must be defined in the median movement.*

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*Please find the exact procedure and the correct technique to attach the facebow and the proper maintenance and care in our AXIOQUICK facebow manual.
This manual can be purchased on request by SAM or can be downloaded as PDF file.*

AXIOCOMP®

Der erste Digitale SAM Artikulator *First Digital SAM Articulator*

vom Erfinder / by inventor of

AXIOGRAPH®, AXIOTRON® - AXIOQUICK® Recorder

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Indispensable in Didactics and Teaching

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Many 3D Views and selectable Features

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Insert the Tilt Support Rods into the holes on the back of the articulator posts.

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They support the articulator in tilted positions in an angle of approx. 40°.

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This improves the view of the anterior teeth.

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Occasionally apply silicone spray, to maintain the O ring (ART 526002).

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Push the cup of the vertical support rod onto the thumb screw of the incisal table.

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It supports the upper part of the articulator when the articulator is open.

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This improves the view of the occlusion, and makes it easier to work on the model (e.g wax up).

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Occasionally apply silicone spray to maintain the support rod and the inner O ring (ART 525002).

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A special feature of the SAM® Plate System is the 15 mm offset arrangement of the mounting plates in the articulator which takes into account the normal lower jaw distal relationship to the upper jaw. By this distal offset of mounting plates position, the mounting of the lower jaw model is easier.

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The different SAM Mounting Plate Systems are available for all SAM articulators and equipment. All SAM instruments can be converted without difficulty.

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The yellow mounting plate is secured with a Thumb Screw and this provides maximum retention for assembly of models with small amounts of plaster. The mounting plates are intended for single use and for storage and easy reattachment of the mounted models.

The low cost of the mounting plates indicates this is the procedure of choice. Not to mention the saving of the costly time to try and recycle the plates. The positioning of dental cast in the articulator is carried out very precisely by means of two precision alignment pins.

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Benefit:

the patient model remains from beginning to end in the same articulator

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The MPS Magnetic Mounting Plate System is an alternative to the historical screw mounting plate. This Magnetic Plate System allows for the quick removal and replacement of the models in the articulator. The conversion of screw type to magnetic system can be done easily and quickly with the appropriate kits (MSP 101K and 101K MPS).

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The AXIOSPLIT® system is used for trouble-free and accurate exchange of models between different articulators (for example, old and new devices), in practice and laboratory. So the models need only to be sent. The articulators stay in place.

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A Master Check ensures the synchronizing with all SAM devices. The mounted models can be interchanged without loss of fitting accuracy in all SAM® devices. You will initially need the AXIOSPLIT Starter Kit with Master Check, norm incisal pin and related matrix holder and an articulator kit for each additional instrument. The exact procedure and more information is described in our AXIOSPLIT-Manual.

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Benefit: Highest precision and your own synchronizing

Screw Type Mounting Plates Attachment

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The Screw Type Mounting Plate has a round hole behind the threaded nut and an elongated slot anteriorly to assure precise plate positioning on all SAM instrumentation. set to the precision mounting plate pins. The elongated hole shows to the front.

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
With one hand the plate is held in position and with the second hand thumb and fore finger the thumb screw is tightened to contact. It is further tightened firmly only 1/4 of a turn. tightened with the knurled screw on the other side.

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The mounting plate is fully seated when it sets firmly on the border areas without gap.

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Take note: The mounting plate is screwed firmly to the articulator, the mounted model may not be set easily and quickly removed on and off as like with the MPS system.

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Note: The mounting plate is disposable and should not be used more than once since the difficulty of removing set gypsum and the cleaning of the mounting plate normally will destroy plate and destroy it's accuracy.

MPS-Magnetic Mounting Plates Attachment

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The adhesion plate is set with the rounded edge downwards into the provided hole.

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
The MPS-mounting plate is attached with the round side containing the slot anteriorly onto the positioning pins.

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Take note: The slightly rounded and shiny side shows downwards when it is inserted into the mounting plate (the rounded edge side shows space for plaster).

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To reproduce or build tooth guidances there should to be an incisal table. It must be able to produce the inclination and direction of the tooth guides.

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
To simulate the anterior teeth and canine guidance authentically, the incisal table is always attached to the upper part of the articulator.

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The incisal table is located as a guiding surface on the top member of the articulator. The flat plane incisal table is used to secure the vertical dimension and to support the incisal pin and can be individualized with plastic via the lateral wings.

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
The incisal pin should be located in the lower part of the SAM articulator.

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
The zero mark on the scale in the middle of the incisal pin is normally the starting position for the dental laboratory.

The incisal pin can be fixed in the restorations bite height. The incisal pin can be raised to 11 mm or lowered up to 6 mm.

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Due to the articulator construction the desired value must be set twice at the incisal pin while raising the occlusion, so e.g. in an increase of 2.5 mm in the occlusion the adjustment at the incisal pin is +5 mm.

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When lowering of the bite to 1 mm, the incisal pin should be lowered 2 mm.

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The micro adjustable incisal pin can be adjusted accurately and within 0.05 mm range. When lateral knurled screw is opened the pin can be finely adjusted on the upper black scale of incisal continuously.

This can be finely adjusted in the micrometer range.

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*To Raise, turn Clockwise.
Rotate Counterclockwise to Lower.*

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Open the lateral knurled screw slightly (about 1/8 turn).

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Rotate the black micro adjustment screw clockwise until the horizontal mark on the scale is in the center of the pin and shows "2" with the vertical mark on the upper black scale "3".

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
The position is fixed with the lateral knurled screw.

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
The axis orbital plane in the patient corresponds in the articulator to the line between the bearing surface of the incisal table to the incisal pin (1) and the centers of the condylar balls (hinge axis) (2).

Cranium-related mounted models in the articulator are an essential part of diagnosis and therapy. The correct position of the maxillary model to the reference plane is an absolute prerequisite for model mounting with a centric registration, but also for setting condylar pathways and even when using mean values.

The Frankfort Horizontal (FH) Plane represents a reference plane for the patient and the articulator. The reference points of the Frankfort horizontal (Porion and orbitale) are identical to the axis-orbital plane (AOE).

*The scissors mechanism of the SAM * anatomic facebow creates and allow for an optimal positioning of the transfer fork (from a functional and aesthetic point of view).*

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Please find the exact procedure and the correct technique to attach the facebow and the proper maintenance and care in our AXIOQUICK facebow manual. This manual can be provided on request by SAM or can be downloaded as PDF file.

www.sam-dental.de

Model mounting in the articulator

The transfer stand AX is used for an easy and fast maxilla mounting in all SAM articulators and laboratory cast mounting instruments. After the transfer the anatomic facebow can remain in the treatment room, only the transferfork assembly is needed for the model transfer. All versions of the SAM transfer fork assemblies are compatible with the transfer stand AX. With the transfer stand the upper jaw position of the patient can be quickly and precisely transferred to the articulator.

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A more detailed manual for model mounting can be provided on request by SAM or can be downloaded as PDF file.

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The occlusal plane indicator set (ART 385) with a plane and a curved plate with Speescher curve can be used for the average-maxillary model mounting or for a prosthetic setting up.

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
The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm.

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The height of the stage device is set to 55 mm (corresponding to the halving of the height - 110 mm between upper and lower part) The medium order setting the depth is 75 mm The angle scale can be adjusted between 0° - 45°, we recommend that the inclination is usually between 0° - 15°.

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is used in the production of full dentures and as an aid for a setting up of the maxillary teeth. It provides a Speesche curve. The marks on the calotte indicate auxiliary lines as the middle of jaw ridge. In the disolved state, the cap may also adjust the inclination angle of the remaining teeth can be determined as occlusion plane.

Cleaning and Maintenance of the Articulator

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Clean the centric lock and the condylar housing with cleaning solvent to eliminate residues of grease, metal and dust.

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The cleaning solvent and the dirt is wiped off with a dry cloth.

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
Lubricate the moving parts with silicone grease (ART 193)

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
Then spread the lubricant while moving in the centric lock.

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*Every now and then, the axis of the condylar housings should be freed from dirt and grease residues, so that the housing remains smooth-running.
Remove successively the thumb screws (A) for the inclination of the condylar housings, the washer (B) and the brass pressure piece (C).*

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The locking screws can be unscrewed with a screwdriver.

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The condylar housing can now be removed by pushing and pulling.

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
All metal parts are cleaned with cleaning solvent and freed from debris.

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After cleaning, all moving metal parts should be lubricated with a silicone grease e.g. like this aluminium shaft.

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While rotating the condylar housing it is inserted into the shaft hole, and snaps into place.

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If the condylar housing does not slip in, press the securing bolt (secured with a spring) down for example with an allen wrench and slide the housing in.

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The brass pressing piece is placed with the beveled side of the condylar axis. After this, insert the washer and insert the thumb screw. Holding this in place with the locking screw.

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All repairs, maintenance, adjustments and engravings or laser inscriptions on your device can be carried out professionally and quickly by our service team.

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For cleaning and maintaining your articulator we recommend our field checking kit (ART 187)

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If you need additional mounting plates or if you have lost a screw, a nut, a washer or any other item from your SAM equipment, please get in contact with us or your local dealer. They have a small stock of most common parts. Don't hesitate to contact them!

You can also download our catalogues for products, spare parts and consumables from our download-page:

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SAM® Präzisionstechnik GmbH ist
*SAM® Präzisionstechnik GmbH is
SAM Präzisionstechnik GmbH