SAM[®] SE SAM[®] SE KIT

HANDBUCH
MANUAL

BIODESIGN - Anatomie als Vorbild

BIODESIGN - Anatomy as benchmark

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
The SAM SE is a patient simulator which allows the end-user to be able to reproduce anatomic patient mandibular movements for examination, diagnosis, and therapy.

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The articulator axis orbital (Frankfort horizontal) plane and all horizontal planes of the articulator are parallel with the table top. Centric relation mounted casts are anatomically related to the patient's cranium with enough space for all types of casts including those with split casts.

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The construction of the SAM SE articulator has the same functional dimensions as all other SAM articulators (Arcon principle).

The upper member guidance surfaces correspond to those found in the cranium as well as the condylar elements and the incisal pin found in the lower part of the articulators relates to the mandible and the lower teeth.

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Anatomically correct based simulation

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
The upper member of the articulator comprised of condylar housings posteriorly and upper incisal guidance table anteriorly. Together, they simulate the TMJ fossa and upper anterior teeth guidance.

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The moving lower member of the articulator with condylar elements posteriorly and the incisal pin anteriorly, together, they simulate the mandibular condylar movement and lower anterior teeth guidance.

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This, therefore, provides a very logical and natural laboratory simulation of the patient movements that is easy to understand and to reproduce.

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The red dot in the center of condylar elements of the articulator lower member correspond to a hinging axis rotation center of the mandibular condyles. Also known as „Hinge Axis“.

The red dotted/dashed line shown above line is the axis-orbital plane (AOP):

The upper anterior incisal table surface along with the posterior condylar element centers of rotation are on the same horizontal plane. The lower incisal pin in contact with the upper incisal guidance table is used to simulate the patient lower anterior teeth incisal edge contact against the lingual surfaces of the upper anterior teeth movement and the resultant normal anterior guidance functional patterns.

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Natural anatomic movement simulation is logical and understandable.

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The SAM articulator is delivered completely assembled and is protected exteriorly with a hard cardboard container which is lined interiorly, on all walls with egg crate foam. (See illustrations).

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When the SAM SE articulator must be transported with or without mounted casts, it is recommended to use ART 598 transport case, optionally available from SAM.

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The transport case is specially designed to receive the SAM SE articulator along with the Axioquick Face Bow and other parts in the SAM SE kits. There is also space for a set of mounted casts.

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
01 articulator lower member (without incisal pin)

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02 articulator upper member (with incisal table)

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05 vertical support rod and 03 incisal pin

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The SAM SE kit comes standard in a transparent plastic case with foam moldings.

This re-packaging for the transport of all parts of the kit has special storage pockets to protect all valuable SAM products during transport.

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- | | |
|--|--|
| 01 articulator lower member (without incisal pin) | 07 individual pockets in the white foam for 3 pairs of mounting plates |
| 02 articulator upper member (incl. incisal table) | 09 facebow without assembly or nasion relator |
| 16 acrylic transfer fork support and 03 incisal pin, | 13 transfer fork assembly / 15 transfer fork |
| 06 Transfer stand SE shown from a front view with the black dove tailed male face bow attachment on the left side, in position for attachment of the transfer fork assembly female counter part. | 11 nasion relator / 05 vertical support rod |
| | X Please leave blank here (antagonistic product) |



Abb. zeigt SE Kit-NT (ART 616MK)
 Fig. shows SE Kit-NT (ART 616MK)

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Contents

✓ enthalten / included
 ✗ nicht enthalten / not included

SE Articulator Kit-AX	SE Articulator Kit-WT	SE Articulator	
✓	✓	✓	01 SAM SE Unterteil / <i>SAM SE Articulator lower member</i>
✓	✓	✓	02 SAM SE Oberteil / <i>SAM SE Articulator upper member</i>
✓	✓	✓	03 Inzisaltift SE - komplett / <i>incisal pin SE, complete</i>
✓	✓	✓	04 Inzisaltisch, plan / <i>incisal table, flat</i>
✓	✓	✓	05 Abstützstift / <i>vertical support rod</i>
✓	✓	✗	06 Transferstand SE / <i>Transfer stand SE</i>
✓	✓	✗	07 8 Stck. MPS-Magnet-Montageplatten / <i>8 pcs. MPS magnetic mounting plates</i>
✗	✗	✓	07 2 Stck. MPS-Magnet-Montageplatten / <i>2 pcs. MPS magnetic mounting plates</i>
✓	✓	✓	08 Haftplatten für MPS / <i>adhesion plates for MPS</i>
✓	✓	✗	09 Transferbogen AXIOQUICK® III / <i>AXIOQUICK® III transfer bow</i>

SE Articulator Kit-AX	SE Articulator Kit-NT	SE Articulator
✓	✓	✗
✓	✓	✗
✓	✗	✗
✗	✓	✗
✗	✓	✗
✓	✓	✗
✓	✓	✗
✓	✓	✗
✗	✗	✓

- 10 Nasionstütze AX - S / *nasion relator AX – S*
- 11 Hygieneschutzkappen / *blue hygienic earpiece caps*
- 12 Bissgabelträger AX / *transfer fork AX*
- 13 Bissgabelträger I-NT / *transfer fork I-NT*
- 14 Bissgabeladapter NT / *transfer fork adapter NT*
- 15 Bissgabel AX / *transfer fork AX*
- 16 Acrylblock / *transfer fork support, acrylic XL*
- 17 Transportkoffer mit Einsätzen / *transparent packing*
- 18 Kartonverpackung / *cardboard packaging*

✓ enthalten / *included*
 ✗ nicht enthalten / *not included*

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
*Horizontal Condylar
Inclination (HCI) scale
in degrees which is
related to the Axis
Orbital Plane.*

Kondylargehäuse
entspricht der
Gelenkpfanne
*condylar housing
assembly to simulate
cranial fossa*

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The incisal pin sleeve is inserted into the lower member until it stops and fixed with the bottom thumb screw.

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We recommend that you attach the incisal pin on the articulator lower member, as this is corresponds to the natural motion simulation.

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The incisal table is inserted into the receiving slot of the articulator upper member.

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
The slots on the articulator upper as lower part are identical in size. This also allows for reverse mounting of the incisal pin and table.

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Preparation - Joining articulator upper with the lower member

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
The centric locking screw on top of the condylar housing can be opened so that the centric lock is movable.

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
The condylar housings are gently lowered onto the condylar elements of both sides.

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With one hand, the lower part of the articulator is held and the upper part pushed back slightly and gently to connect both parts.

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
Pressure is applied on the condylar housing. You will hear a clicking noise as the centric lock snaps over the condyle balls, the upper member is now connected to the lower.

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
The pressure spring automatically keeps the upper member of articulator in centric position.

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
The centric lock is closed when the spring is fully extended as shown. The position of the bolt can be fixed with the centric lock screw.

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Caution: Do not touch the end of the spring bolt with your thumb!

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The centric lock is open when the spring is compressed. The position of the bolt can be fixed with the centric lock screw.

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All SAM® articulators are available either with a thumb screw or with magnetic mounting plates system MPS. For the attachment of the AXIOSPLIT® system the screw thread version is needed. A conversion from one cast mounting system to the other via different additional kits is easy to handle for users (See page 23).

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MPS MAGNETIC MOUNTING PLATES YELLOW

- schneller Modellwechsel im Artikulator
quick attachment and removal of the cast
- sicheres Reponieren der Modelle
precise positioning of the casts

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
Screw type mounting plates attachment

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The Screw Type Mounting Plate has a round hole behind the threaded nut and an elongated slot anteriorly to assure precise plate positioning on all SAM instrumentation. set to the bolts, the rounding with the elongated hole shows to the front.

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Take note: The mounting plate is screwed firmly to the articulator, the model may not be set quickly in and out like with MPS system.

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With one hand the plate is secured, with the second tightened with the knurled screw on the other side.

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
The mounting plate is fully seated when it bears on the border areas without gap.

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Note: *The mounting plate is disposable and should not be used more than once for reasons of quality.*

MPS-Magnet Montageplatten-Befestigung

MPS-magnet mounting plates attachment

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The adhesion plate is set with rounded edge into the provided hole.

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Take note: The slightly rounded and shiny side shows downwards when it is inserted into the mounting plate (the edgy side shows page for plaster).

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The MPS-mounting plate is attached with the round side forward into the positioning pins.

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Note: *The mounting plate is an disposable and should not be used more than once for reasons of quality.*

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To facilitate the MPS mounting plate removal from the upper member a plate ejector is installed.

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Pressing the eject button („push“) the mounting plate is slightly raised in the front and can be easily removed.

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To remove the mounting plate from the bottom, the incisal pin is useful.

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Here is pressed gently from the bottom rear of the mounting plate. Caution - Do not use sharp pointed objects!

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Converting Screw type plate to MPS-System

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Use appropriate Allen Wrench to loosen set screw holding mounting plate Thumb Screw in position.

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Position magnetic cup threaded screw to position as shown in illustration above.

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Use yellow combination wrench to loosen MPS magnet cup nut that is holding MPS cup in place. Save all parts for later use.

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Place threaded screw type piece into position as shown. Make sure to put plastic washer in place. Align Thumb screw to position as shown and place on top of threaded extension. Use Allen wrench to tighten set screw firmly

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More detailed instructions are in the respective conversion kits
(Order-Nr.: ART 114K and MPS 101K).

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The design of the SAM SE can be tilted backward and placed in an inclined position.

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
The SAM SE - is like all SAM articulators are designed - originally in an inverted position as well as upright.

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Notice how this improves the look at the occlusion while observing and controlling the function and effect of the condylar elements of the lower articulator member.

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Simply place the cupped shaped part of the vertical support rod onto the upper thumb screw.

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The vertical support rod holds the articulator in a stable horizontal position when open.

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Notice how this improves the view of the occlusion

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You maintain the support rod best in which the inner O-Ring (Art 525002) occasionally lubricated with silicone spray.

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
*often called as **Bennett angle** is a Mediotrusion curve and no straight pathway, as in most mean value articulators. The Bennett guide is already given in the condylar housing and corresponds to the most common patient situation in a healthy TMJ.*

It can not be set individually.

In a full adjustable SAM articulator (e.g. SAM 3) this pathway corresponds to the green Bennett use at 5° individual setting.

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
*Learn more about Bennett and lateral movements in our end user training
(information at www.cde-info.de)*

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





When the mandible makes a three-dimensional movement then the SAM® AXIOGRAPH® or AXIOQUICK® Recorder will show a curved mediotrusion pathway in both frontal and horizontal planes. This curve corresponds to the most common of all TMJ movement and is used to adjust the articulator.

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
The „successor“ to the AXIOGRAPH® is the AXIOQUICK® Recorder. It is a non-contact measurement system using ultrasonic technology.

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The Bennett movement is a downward, inward, forward movement or downwards, forward, inward movement all at the same time of the condyle on the mediotrusion (balancing, non-working) side while the laterotrusion (working, rotating) side moves virtually straight lateral for about 0.3 mm.

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An exaggerated mediotrusion movement is shown here. The upper and lower jaws move in the opposite direction. The canine guidance causes disclusion of the posterior teeth when the mandible moves laterally and/or in a protrusive direction.

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To perform a laterotrusion movement, the locking screw of the spring-loaded bolt is loosened.

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By holding the condylar element firmly against the medial wall of the condylar housing and pushing the condylar element downward, inward, and forward the condylar element guiding on the condylar housing pathway will move the lower member towards the midline to perform the mediotrusion movement.

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
Normal opening/closing patient condylar movement without translation is a simple rotation of the condyle head within the disc held firmly against the eminentia. Maximum opening/closely, beyond a simple rotation of the condyle, is an immediate downward (Z axis) movement followed by a forward (X axis) translation which causes the disc to rotate backwards on the head of the condyle and forwards under the eminentia. This is simulated in the articulator with the condylar pathway inclination and the joint path curvature.

The individual values of the mandible movement can be recorded.

Once you have the Axiographic values it is a simple procedure to select the appropriate condylar housing pathway inserts for Bennett

Guidance or Condylar pathway inclination and curvatures but you must have a SAM 3 articulator. In SAM SE articulator the condylar inclination (angle) can be adjusted, but not the joint path curvature. For the joint path curvature is a statistical average applied after anatomical ideal. SAM SE articulator has equivalent of the Condylar insert 1 (white, low curvature along with the Green Bennett guidance pathway set at maximum Bennett angle of 10 degrees.)

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

The simulation of the patient pathway is a combination of three parameters:

- *Condylar ball diameter*
- *Condylar pathway curvature*
- *Condylar pathway inclination*

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The path shown is the real path of the center of rotation of the temporomandibular joint head of the patient.

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To simulate this as faithfully as possible, the condyle of the envelope of condylar pathway elements (joint path curvature) must follow.

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The joint path is most often a curved trajectory. This was reflected in the condylar SAM SE.

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In addition the condylar inclination must be adjusted.

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Setting the condylar pathway inclination

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For this the condylar housing is loosened or tightened with the rear thumb screw (arrow).

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

The condylar inclination can now be set at 15° - 65° with respect to the axis-orbital plane.

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In addition to the set horizontal condylar pathway inclination (HCI) is given following value by the joint path curvature:

- 1 mm: + 6,0°
- 2 mm: + 6,0°
- 3 mm: + 5,0°
- 5 mm: + 4,0°
- 7 mm: + 2,5°
- 10 mm: + 0,0°

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Simulation on the protrusive movement

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
To simulate a protrusive movement, the locking screw of the spring bolt (arrow) is loosened and the condylar housing shifted slightly to the rear.

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This achieves a flowing movement that moves on both sides as protrusion. The mandible can move forward.

Important: Do not block the spring pin with your thumb.

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With protrusion inserts (optional) in 1 mm steps (1-6 available), a protrusive position of the lower member in relation to the upper member can be fixed anteriorly.

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These can be used in all SAM SE articulators (from April 2015)

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

Cranium-related mounted models in the articulator are essential part of diagnosis and therapy. The correct position of the maxillary model to the reference plane absolute prerequisite for model mounting with a centric registration, but also for setting condylar pathways and even when using mean values.

The Frankfort Horizontal (FH) represents a reference plane for the patient and the articulator. The reference points of the Frankfurt horizontal (Porion and orbitale) are identical to the axis-orbital plane (AOE).

The scissors mechanism of the SAM[®] anatomic facebow creates and allow for an optimal positioning of the transfer fork (from a functional and aesthetic point of view).

The axis orbital in the patient corresponds in the articulator to the line between the bearing surface of the incisal table to the incisal pin (1) and the centers of the condyle balls (hinge axis) (2).

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
Hard silicone (eg metal-Bite of R-Dental) is applied so that the bore holes on the bottom are masked and then the top is coated evenly thin.

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With minimal pressure the coated transfer fork is pressed on the teeth parallel to the occlusion, without pressing through the material. The bite fork middle and the anterior center should coincide.

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After hardening the transfer fork is removed and all the deep impressions, shaped fissures and interdental spaces were removed. Only the impressions of cusp tips are to be shown!

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
The transfer fork is inserted again and supported with cotton rolls in the area of the premolars. The patient closes the lower jaw and the mandibular teeth support the transfer fork.

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The blue hygienic earpiece caps are placed on the black earpieces. Placed correctly the notch shows outside, the flat side towards the inside.

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The nasion relator (ATB 394) is positioned correctly on the facebow (ATB 303), when it is fully inserted between the pins on the bottom.

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The transfer fork assembly is pushed onto the attachment at the bottom of AXIOQUICK® facebow to positioning pin. The thumbscrew is tightened by half a turn and is visible to the dentist.

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The attachment connection is secured by tightening the knurled screw. The other clamp-connections are to be loosened only a half turn.

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The facebow is held by patient with the thumb and index finger on the angle to the ear pieces.

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The patient inserts the earpieces with hygienic earpiece cap into the ear canal to slightly block the hearing. A humming sound is made to equalize the earpiece position in the ears. The earpieces are then moved straight upwards and inwards to contact the top inner wall. And then the hearing is fully restored.

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Notice initial positioning in the ear.

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To bring the face bow in the exact position, the patient raises the earpieces straight upwards and then inwards into a small superior depression onto anatomic porion. The patient can hear clearly again.

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
With earpieces held firmly in the proper position, the nasion relator is brought into light contact with the nasion and fixed in position without pressure using the thumb screw.

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All clamps and the thumb screws are loosened. The transfer fork rod is seated firmly to position and the thumb screw is tightened.

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The interpupillary line levelling rod is used to control the parallelism of the interpupillary line and facebow. The facebow is now esthetically correct and aligned to the Frankfort Horizontal and the interpupillary line.

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With one hand, the transfer fork assembly is stabilized, with the other hand, the screw of the NT-clamp is tightened.

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***An incorrectly applied
facebow leads
inevitably to an
incorrect transfer.***

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The toggle on the horizontal rod is tightened after all other adjustments have been made and set.

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
To remove the facebow the patient opens his mouth, while the user presses on the front of the bow.

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The patient is instructed to open the mouth and the face bow assembly is removed from the patient.

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
*The red arrows mark the rods and holes of the transfer fork which must be absolutely lubricant and dirt free for. conservation of the friction.
Degreasing with alcohol assures a better clamping action.*

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The blue arrows mark internal and external threads of the transfer fork assembly which should periodically lubricated with dental hand piece oils or similar types of lubricants.

Reinigung, Desinfektion und Sterilisation *Cleaning, disinfection and sterilization*

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Facebow and transfer fork assembly is designated as a „non-critical medical device“, so a wipe disinfection is recommended.

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Facebow, transfer fork assembly, adaptors and nasion relator should NEVER be placed in chemical solutions or ultrasound.

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Cleaning, disinfection and sterilization

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*The blue hygiene protection caps are designed **only for the single use.***

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The transfer fork ATB 396 should be sterilized in the Autoclave

The transfer fork is classified as a „semicritical medicine product“ (according to MPG/act concerning medicine products) and must be sterilized as such. For this the steam sterilization is recommended in the fractionated vacuum procedure for 134 degrees centigrade for 10 minutes in a device after German Institute for Standardization DIN EN ISO 17665-1.

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Remove the incisal pin and replace it with the transfer stand. The Centric lock must be closed.

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The transfer fork is slid onto the male adapter all the way onto the dovetail attachment of the transfer stand and secured with the thumb screw.

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The transfer fork should not change its position. Therefore, it should always be supported.

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The SAM acrylic block can be used in conjunction with silicone or plaster to support the transfer fork.

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Estimate the required amount of plaster by closing upper member of the articulator.

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The creamy mixed mounting plaster is first applied to the upper mounting plate filling all retention areas of the plate. Gypsum should flow freely.

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A small stacked amount of creamy gypsum is placed on top of a previously moistened upper cast to assure adhesion.

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Carefully close the upper member onto the stacked opposing gypsum.

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If the lower cast is mounted with the articulator locked in Centric then it is recommended to remove the lower incisal pin when closing the articulator together to prevent gypsum contact with the incisal pin.

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The setting on the incisal pin is obtained from the measurement result of the thickness of the centric Registration in the premolar region multiplied by two.

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With a centric registration the lower jaw model is placed onto the maxillary model.

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To prevent a change in position during the assembly of models, the two models should be fixed with hotmelt or wax.

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The acrylic block compensates for the angle of the mandible. Thus, the model surface is oriented parallel to the table plane.

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Creamy gypsum is placed onto the lower mounting plate in the same manner as was used for the upper plate. The model is slightly dampened, so that the fresh plaster blends well with the model.

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The articulator is closed and is not moved until gypsum has set.

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The incisal pin and incisal table adjusted back into position.

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The SAM SE is an easy to care for articulator. Plaster residues can be removed with water, soap and a soft brush.

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To remove greasy stains, we recommend mild grease solvents.. Condylar ball and housing should always be clean and kept dry (oil and grease free).

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The small centric screw on condylar housing, and the set screws on the incisal pin and table should be occasionally cleaned with mild grease solvents. and then greased (eg with silicone spray or Molykote).

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It is important to ensure that the centering pins for mounting plate are clean to ensure an exact fit. Also, all mounting plate contact surfaces should be free of gypsum and other debris.

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Avoid air drying liquid gypsum onto the articulator. This can damage the instrument and make it unrepairable.

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Prior to use, make a quick check of the SAM SE by viewing the action of the lower incisal pin into the slot of the upper member with the incisal table removed.

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ATB 338

Acrylblock
transfer fork support, acrylic XL

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SVst – perfect for entry into the world of loupe glasses

Reasonably priced, lightweight and sporty loupe glasses for students. You will discover a new world of vision with 2.5x magnification and a very wide field of vision.

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SV – the perfect synthesis of form and function

The quality of these loupe glasses is reflected in a highly efficient design and minimal weight. The optics ensure a very wide field of vision and outstanding depth of field at a 2.7x magnification. The ergonomic and practice-oriented design mean these loupe glasses are a tool that will make your work easier.

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The tried and tested optics offer a bright and brilliant image with no chromatic aberrations. Great depth of field and excellent imaging quality are impressive features in challenging situations and make this versatile magnifying glass an indispensable tool for one or more users in everyday practical situations.

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The compact and portable starLight nano LED lighting system optimally lights up your work area. The starLight nano is easy to use and can be adjusted for use with starVision loupe glasses or any other standard loupe glasses.

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SAM® Präzisionstechnik GmbH